



Key Code TSMX7779B  
www.oxid.com/ifu

Europe + 800 135 79 135  
CA 1 855 805 8539

US 1 855 236 0910  
ROW +31 20 794 7071

# PathoDx Strep Grouping Extraction Reagents


<b>REF</b>	R62050.....Reagent 1 (7 ml)
	R62052.....Reagent 1 (14 ml)
	R62055.....Reagent 2 (7 ml)
	R62057.....Reagent 2 (14 ml)
	R62060.....Reagent 3 (14 ml)

## INTENDED USE

Remel PathoDx™ supplementary Strep Grouping reagents are intended for use in the PathoDx Strep Kits, which contain all components necessary for the various systems. These reagents are provided as replacement items. **For complete information and correct use of this product, please refer to the Instructions For Use which accompanies the PathoDx Strep Grouping kits (R62025, R62076 and R62150).**

## DESCRIPTION, PREPARATION FOR USE AND RECOMMENDED STORAGE CONDITIONS

See also **Warnings and Precautions**.



Store at 2 to 28°C. If crystallization of Extraction Reagent 3 should occur, allow to reach room temperature before use.

### EXTRACTION REAGENT 1

#### Reagent 1 (Red cap)

One dropper bottle containing 7.0 ml (R62050) or 14 ml (R62052) of extraction reagent. Store tightly capped at 2 to 28°C. Use on or before the expiration date marked on the label.

### EXTRACTION REAGENT 2

#### Reagent 2 (Blue cap)

One dropper bottle containing 7.0 ml (R62055) or 14 ml (R62057) of extraction reagent. Store tightly capped at 2 to 28°C. Use on or before the expiration date marked on the label.

### EXTRACTION REAGENT 3

#### Reagent 3 (Green cap)

One dropper bottle containing 14 ml (R62060) of neutralization reagent. Store tightly capped at 2 to 28°C. Use on or before the expiration date marked on the label.

## WARNINGS AND PRECAUTIONS



The reagents are for *in vitro* diagnostic use only.

For professional use only.

Please refer to the Safety Data Sheet and the product labeling for information on potentially hazardous components.

## HEALTH AND SAFETY INFORMATION

- Extraction Reagent 2 contains Acetic Acid which is classified per applicable European Economic Community (EEC) Regulations as corrosive. The following are the appropriate Hazard (H) and Precautionary (P) statements.
- The mixture of Reagents 1 and 2 has an acidic pH and should be considered potentially hazardous until neutralized with Reagent 3. While these reagents are not likely to be harmful to intact normal skin, special care should be taken to avoid contact with eyes, sensitive mucous membranes, cuts and abrasions. If skin contact does occur, wash the site with soap and water. In case of eye contact, rinse liberally with water.

### EXTRACTION REAGENT 1

## WARNING




H302	Harmful if swallowed.
H319	Causes serious eye irritation.
P264	Wash face, hands and any exposed skin thoroughly after handling
P301 + P312	IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### EXTRACTION REAGENT 2








## DANGER



H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves/ protective clothing/eye protection/face protection.
P301 + P330 + P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P303 + P361 + P353	IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/shower.
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

EXTRACTION REAGENT	3	H319	Causes serious eye irritation.
<b>WARNING</b> 		P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

**SYMBOL LEGEND**

	Catalogue Number
	<i>In Vitro</i> Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Manufactured by



IFU X7779B, Revised March 2016      Printed in the UK

 Remel Europe Ltd. Clipper Boulevard West,  
 Crossways, Dartford, Kent, DA2 6PT UK

For technical assistance please contact your local distributor.