

COX2

Catalog # RM-9121-S0, -S1, or -S (0.1ml, 0.5ml, or 1.0ml)

Catalog # RM-9121-R7 (7.0ml)

lease note the data sheet has been changed effective Apr 17, 2009

INTENDED USE:

- **For In Vitro Diagnostic Use:** This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.
- **Description:** COX2 (Cyclooxygenase-2) is an inducible enzyme. It is involved in the response of cells to growth factors, tumor promoters, and cytokines that induce its expression. Given its role in synthesizing prostaglandins, COX2 is therefore of interest in studying immune response regulation. COX2 is induced by a wide variety of stimuli and was initially identified as immediate-early growth response gene. In addition, COX2 expression markedly increased in 85-90% of human colorectal adenocarcinoma whereas COX1 levels remain unchanged.
- **Expected Staining Pattern:** Cytoplasmic and membrane
- **Positive Control:** Lung and colon CAs

MATERIALS PROVIDED:

COX2 (refer to catalog number):

- #RM-9121-S (or -S0, -S1): Tissue culture supernatant, concentrated, with 0.09% Sodium Azide.
or
- #RM-9121-R7: (7.0ml) of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing stabilizing protein and 0.015mol/L sodium azide.
- **Antibody Concentration:** Not known
- **Host:** Rabbit
- **Mol. Wt. of Antigen:** 70kDa
- **Epitope:** C-terminus
- **Species Reactivity:** Human, Mouse, Rat. Others not tested.
- **Clone Designation:** SP21
- **Ig Isotype / Light Chain:** Rabbit IgG
- **Immunogen:** A synthetic peptide from C-terminus of rat cox2
- **Microbiological State:** This product is not sterile.

MATERIALS REQUIRED, BUT NOT PROVIDED:

- **Antibody Diluent:** For concentrated antibodies, the antibody must be diluted before using. Use Lab Vision Antibody Diluent (catalog # TA-125-UD). Refer to diluent product instructions for use.
- **Negative Control Reagent:** Refer to the "General Protocol" instructions.
- **Visualization System:** Refer to the "General Protocol" instructions.

METHODS AND PROCEDURES:

Specimen Preparation	Refer to the "General Protocol" instructions.
Dilution of Concentrated Antibody	1:50-100 in antibody diluent using UltraVision LP Systems
Tissue Section Pretreatment	Staining of formalin-fixed tissue sections requires treating the tissue sections in boiling 10mM citrate buffer, pH 6.0 (Lab Vision catalog # AP-9003), for 10-20 minutes followed by cooling at room temperature for 20 min.
Primary Antibody Incubation Time	20 minutes at Room Temperature
Visualization	To detect antibody, follow the instructions provided with the visualization system.

STORAGE and STABILITY:

This product contains sodium azide and is stable for 24 months when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent is not stored as recommended, performance must be validated by the user.