

## Rheometer accessory

# Rheological measurements of UV curing reactions at ambient temperature

## Discover the cost-efficient entry model

**Author**

Philipp Beutler and Cornelia  
Küchenmeister-Lehrheuer  
Thermo Fisher Scientific,  
Karlsruhe, Germany

**Keywords**

HAAKE MARS 40/60 and HAAKE  
MARS iQ Air Rheometer, measuring  
plate cover, UV light source

For measurements at ambient temperature under defined exposure to UV light a simple measuring cell is available which can be used in combination with a Thermo Scientific™ HAAKE™ Rheometer. It is characterized by a simple design and is suitable as a very cost-effective model for beginners in the field of UV measurements. Measuring cells for UV assisted thermal curing<sup>1</sup> or UV measurements with simultaneous structure analysis using FT/IR spectroscopy<sup>2</sup> extend the accessory portfolio of Thermo Scientific™ HAAKE™ HAAKE™ MARS™ Rheometers.

This CM-OP-P module (for Combined Module - OPTical -Plate) consists of a “triangular” plate, on which a holder for a glass plate and a port for a UV light collimator below that, are mounted (Figure 1). The CM-OP-P module does not offer temperature control, but is equipped with a temperature sensor.

A commercially available UV light source<sup>3</sup> can be used, and controlled via the Thermo Scientific™ HAAKE™ RheoWin™ Rheometer software (Figure 2).

Plate/plate measuring geometries with diameters up to 20 mm and from various materials (e.g. titanium, or aluminum as disposable version) are available.



Figure 1. CM-OP-P module for UV measurements at ambient temperature.




Figure 2: Curing curves as function of different UV light intensities.

## Ordering information

Product	Cat. No.
CM-OP-P module for cost effective UV measurements at ambient conditions	222-2498
Rotor P20/Ti with "Connect Assist" and ceramic shaft	222-2090
Adapter P3 for disposable plates D PXX /Al with "Connect Assist" and ceramic shaft	222-2290
Disposable geometry D P20/Al (40 pcs.) made of aluminum	222-2154
Set of glass plates for UV measuring cell (5 pcs.)	222-2388

## References

1. F. Meyer, Universal Peltier temperature module for UV-curing with mercury vapour lamp, Thermo Fisher Scientific Product information P064
2. C. Küchenmeister-Lehrheuer, J. Nijman, Measuring cell for UV assisted thermal curing and simultaneous FTIR spectroscopy, Thermo Fisher Scientific Product information P037
3. Ph. Beutler, Mercury-vapor lamp or LED?, Thermo Fisher Scientific Product information P072

Learn more at [thermofisher.com/rheometers](https://thermofisher.com/rheometers)

thermo scientific

For research use only. Not for use in diagnostic procedures. For current certifications, visit [thermofisher.com/certifications](https://thermofisher.com/certifications)

© 2025 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. P088 09/25