

HPLC & UHPLC

List of compendial methods

Thermo Scientific Charged Aerosol Detectors

Charged aerosol detection is a reliable technology that will change the way you view every sample. The Charged Aerosol Detector (CAD) can detect all non-volatile, and many semi-volatile analytes, with uniform response. Charged aerosol detection can be used for the analysis of pharmaceuticals (large and small molecule), biomolecules, foods and beverages, specialty chemicals, and polymers.

In the following tables are examples of standardized methods—for example, American Society for Testing and Materials (ASTM International), Chinese Pharmacopoeia (CP), European Pharmacopoeia (EP), International Organization for Standardization (ISO), and United States Pharmacopoeia (USP)—using charged aerosol detection technology.

General chapter information on Charged Aerosol Detectors

Chinese Pharmacopoeia

"The latest HPLC developments and application progress are fully detailed. Information about multidimensional HPLC, charged aerosol detection, adjusted chromatographic conditions, and common qualitative analysis methods was added." 2020 edition.

European Pharmacopoeia

"Ultraviolet/visible (UV/Vis) spectrophotometers (including diode array detectors) (2.2.25), are the most commonly employed detectors. Fluorescence spectrophotometers, differential refractometers (RI), electrochemical detectors (ECD), light scattering detectors, charged aerosol detectors (CAD), mass spectrometers (MS) (2.2.43), radioactivity detectors, multi-angle light scattering (MALS) detectors or other detectors may be used."



Specific monographs using Charged Aerosol Detector

Organ-ization	Analytes	Method Reference	Matrix	Column	Mobile Phase	Gradient or Isocratic	Status
							Additional Information
ASTM	Lipids	ASTM E3297-21 Standard Test Method for Lipid Quantitation in Liposomal Formulations	Liposomal Formulation	RP C18	Aqueous acetonitrile/ methanol	Gradient	Active
							Additional info: AppsLab Library
EP	API and counterions	Validation and application of an HPLC-CAD-TOF/ MS method for identification and quantification of pharmaceutical counterions. Pharmeuropa bio & scientific notes 2014, 81-91	None	Mixed mode Thermo Scientific™ Acclaim™ Trinity P1, 3.0 × 50 mm, 3 µm (P/N 071388)	A) 0.2 M ammonium formate in water (pH 4.0) B) water C) acetonitrile	Gradient	Published for consideration
							Additional info: AppsLab Library Mixed-mode HPLC columns
EP	Aspartic acid	Aspartic acid Pharmeuropa 22.4, October 2010	None	RP C18, 4.6 × 150 mm, 5 µm	4% 0.36 g/L perfluoroheptanoic acid (PFHPA) in methanol 96% 0.36 g/L (PFHPA) in water (%/% v/v)	Isocratic	In process/under consideration
							Additional info: AppsLab Library
EP	Gadobutrol monohydrate	Gadobutrol Monohydrate–Ph. Eur. 9.0 07/2016:1215	None	RP end-capped phenylhexylsilyl silica gel, 4.6 × 250 mm, 3 µm	A) 0.5% acetonitrile and 99.5% water pH 3.6 adjusted with formic acid B) acetonitrile	Gradient	In official text
EP	Ibandronate sodium monohydrate	Ibandronate Pharmeuropa 22.4, October 2010	None	RP octadecylsilyl silica gel with embedded strong anion exchange groups, 4.6 × 150 mm, 2.7 µm	A) water B) 150 mL acetonitrile plus 150 mL of 11.4 g/L trifluoroacetic acid, filled to 1 L with water	Gradient	In process/under consideration
EP	Topiramate	Topiramate Pharmeuropa 27.4, October 2015	None	RP solid core pentafluorophenylsilyl silica gel, 4.6 × 100 mm, 2.6 µm	A) 1.93 g/L ammonium acetate, pH 3.5 B) acetonitrile	Gradient	Adopted for next version
EP	Valine and impurities	Determination of the purity of valine by isocratic liquid chromatography coupled with charged aerosol detection Pharmeuropa bio & scientific notes. 2015; 2015, 11-18	None	RP C18, 4.0 × 150 mm, 3 µm	20 mM Perfluorobutyric acid in acetonitrile/ water (10:90 v/v)	Isocratic	Published for consideration
							Additional info: AppsLab Library
EP	Vigabatrin	Vigabatrin Pharmeuropa 30.2, April 2018.	None	RP end-capped solid core phenylhexylsilyl silica gel, 4.6 × 100 mm, 2.7 µm	19.5% methanol, 80.5% water (v/v) with 2.1 g/L perfluoroheptanoic acid (PFHPA)	Isocratic	Adopted for next version

Organ-ization	Analytes	Method Reference	Matrix	Column	Mobile Phase	Gradient or Isocratic	Status
							Additional Information
ISO	PEG with molecular mass greater than 400 g/mol	ISO 16560:2015(en) Surface active agents—Determination of polyethylene glycol content in nonionic ethoxylated surfactants—HPLC method	Nonionic ethoxylated surfactants that are soluble in methanol or methanol/ water/H ₂ O and have [PEG]>0.1%	Reversed phase (RP) C18, 4.6 × 250 mm, 5 µm	A) water B) methanol	Gradient	In official text
							Additional info: Surfactants Application Notebook
ISO	Triton X-100, octylphenol ethoxylates (CAS 9002-93-1) IGEPA CO-630, nonylphenol ethoxylates (CAS 68412-54-4)	ISO 18254-2:2018 Method for the detection and determination of alkylphenol ethoxylates (APEO)—Part 2: Method using NPLC	Textiles	Normal phase (NP) applying two columns: 1) C18 column, 4.6 × 50 mm, 1.7 µm 2) hydrogen bond adsorption, 4.6 × 150 mm, 3 µm	A) acetonitrile with 0.1% formic acid B) methanol with 0.1% formic acid and 0.01% ammonium formate	Gradient	In official text
							Additional info: AppsLab Library
USP	Deoxycholic acid powder	Deoxycholic (desoxycholic) acid—since USP 40 NF 35 S1	None	Type L1, Thermo Scientific™ Acclaim™ 120 C18, 4.6 × 150 mm, 3 µm (P/N 059133)	A) 0.1% formic acid in water B) 0.1% formic acid in acetonitrile	Gradient	In official text
							Additional info: AppsLab Library Application Note Poster
USP	Metoprolol succinate powder	Content of metoprolol related compound H and metoprolol related compound I, USP 41(3) In-Process Revision: Metoprolol succinate. Proposed change to United States Pharmacopeia and National Formulary USP 38-NF33;	None	Hydrophilic interaction liquid chromatography (HILIC) solid core silica gel with five hydroxyl bonded ligands, 4.6 × 150 mm, 5 µm	85% acetonitrile 15% 0.1 M ammonium formate in water, pH 3.2	Isocratic	In process/under consideration
							Additional info: AppsLab Library Application Note Poster

For more references download the Charged Aerosol Detection bibliography highlighting the breadth and scope of different analytical methods found in the literature, click [here](#).

 Learn more at thermofisher.com/CAD

For Research Use Only. Not for use in diagnostic procedures. © 2019-2022 Thermo Fisher Scientific Inc. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. **FL73052-EN 0922M**

thermo scientific