



Evaluating Technical Suitability of High Resolution Mass Spectrometry for Simultaneous Screening and Quantification of 47 Drugs in Human Urine

Ana Celia Grenier, PhD

- ▶ Introduction
- ▶ Materials and Methods
- ▶ Results
- ▶ Conclusions
- ▶ Questions

Introduction

Why did we consider High Resolution Accurate Mass(HRAM) instrumentation as an alternative to immunoassay?

- ▶ Most drug screening is performed in two steps comprising an initial immunoassay screen followed by confirmation by liquid chromatography mass spectrometry.
- ▶ Monitoring various drugs of abuse using HRAM LC-MS allows for simpler sample preparation, and increased analytical specificity and selectivity.
- ▶ Decrease sample analysis cost while increasing accuracy and the number of analytes screened and quantified on a single injection using HRAM LC-MS.

Introduction

Types of Screening	
EIA/Elisa	LC MS Accurate Mass
<ul style="list-style-type: none"> •Fast Cycle Time (seconds) •Low Sensitivity •Low Specificity •Pos/Neg Results Only – drug families, no individual drugs 	<ul style="list-style-type: none"> •High Sensitivity •Very High Specificity •Screening and Quantitation are Accurate •Cycle time in minutes (not seconds)



Introduction

Types of Screening using Accurate Mass

Full MS Scan	Full MS Scan/MS ²
<ul style="list-style-type: none">•Retention Time•Accurate Mass of the Parent•Isotopic Pattern	<ul style="list-style-type: none">•Retention Time•Accurate Mass of the Parent•Isotopic Abundance Pattern•Accurate Mass of Fragments•Library Search



Introduction

Types of Screening Applications Using Full MS Scan/MS²

Method Screens for hundreds of compounds	Method Screens for a restricted number of compounds
<ul style="list-style-type: none">•Ability to set up screening for hundreds of targeted and non-target compounds•Refer back to the data anytime without reruns to investigate samples for both target and non-target compounds•Calibrators and QCs do not have all the compounds present at once	<ul style="list-style-type: none">•Inclusion list has only targeted compounds to perform MS²•Calibrators and QCs have all the compounds present at once•MS² is collected for all the compounds present in calibrators and QCs•All calibrators and QCs have to pass screening criteria.



Introduction

Reasons for choosing a targeted Screening method

- ▶ Require to have a certified reference standard for screening
- ▶ Demonstrate linearity
- ▶ QCs and calibrators need to be included in the analytical run
- ▶ Results are not reported below LLOQ
- ▶ Screening and Quantitation can be done with one single injection – Calibrators and QCs are required for both



Materials and Methods

Calibrators and quality control (QC) samples were prepared in Surine™ using certified reference standards from Cerilliant™

Sample Preparation Procedure

75µL of sample were diluted with 300 µL a master mix solution containing IMCSzyme®, and 47 internal standards, one for each of the analytes. Samples were incubated at 65°C for 60 minutes, and centrifuged at 4000 rpm for 7 min before analysis.

Materials and Methods

LC Parameters

Instrument: UPLC

Analytical Column: Agilent Poroshell™ 120 SB-C18 2.1 x 50mm 2.7µm

Gradient: A: 0.1 % Formic Acid acid in LCMS grade water
B: 0.1 % Formic Acid in Methanol

Time (min)	A (%)	B(%)
0.50	95	5
1.50	75	25
4.60	40	60
4.61	5	95
5.00	5	95
5.01	95	5

Flow Rate: 0.50 mL/min

End time: 5.5 min

Column Temp: 30 °C

Injection volume: 7.5 µL

Autosampler Temp: 5 °C



Materials and Methods

Mass Spectrometer

Source Parameters

Mode: HESI+

Sheat gas flow rate: 65

Aux gas flow rate: 20

Sweep gas flow rate: 2

Spray voltage: 3.5 kV

Capillary Temp: 300 °C

Aux gas heater temp: 420 °C

S-lens RF Level: 55

Full MS

Scan Range: 130–480 m/z

Resolution: 17 500, 35 000, 70 000

MS/MS

First fixed mass: 50 m/z

Resolution: 17 500

ddMS²

Top 2, 3, 4, 5 (targeted)

First fixed mass: 50 m/z

Resolution: 17 500



Results

ANALYTE	Internal Standard	Analyte Calibration Range	ANALYTE	Internal Standard	Analyte Calibration Range
2-Hydroxyethylflurazepam	2-Hydroxyethylflurazepam-D4	20-1000	Methadone	Methadone-D3	20-1000
6-Acetylmorphine	6-Acetylmorphine-D3	20-1000	Methamphetamine	Methamphetamine-D5	100-5000
7-Aminoclonazepam	7-Aminoclonazepam-D4	20-1000	Morphine	Morphine-D3	20-2000
α -Hydroxyalprazolam	α -Hydroxyalprazolam-D5	20-1000	Norbuprenorphine	Norbuprenorphine-D3	20-1000
α -Hydroxymidazolam	α -Hydroxymidazolam-D4	20-1000	Norcodeine	Norcodeine-D3	20-1000
α -Hydroxytriazolam	α -Hydroxytriazolam-D4	20-1000	Nordiazepam	Nordiazepam-D5	20-1000
Alprazolam	Alprazolam-D5	20-1000	Norfentanyl	Norfentanyl-D5	2-100
Amphetamine	Amphetamine-D5	100-5000	Norhydrocodone	Norhydrocodone-D3	20-1000
Benzoylcegonine	Benzoylcegonine-D3	20-1000	Normeperidine	Normeperidine-D4	20-1000
Buprenorphine	Buprenorphine-D4	20-1000	Noroxycodone	Noroxycodone-D3	20-1000
Carisoprodol	Carisoprodol-D7	20-1000	Norpropoxyphene	Norpropoxyphene-D5	40-2000
Cocaethylene	Cocaethylene-D3	20-1000	O-Desmethyltramadol	O-Desmethyltramadol-D6	20-1000
Codeine	Codeine-D6	20-2000	Oxazepam	Oxazepam-D5	20-1000
EDDP	EDDP-D3	20-1000	Oxycodone	Oxycodone-D3	20-2000
Fentanyl	Fentanyl-D5	2-100	Oxymorphone	Oxymorphone-D3	20-2000
Gabapentin	Gabapentin-D10	20-1000	PCP	PCP-D5	20-1000
Hydrocodone	Hydrocodone-D3	20-2000	Pentazocine	Pentazocine-13C3	20-1000
Hydromorphone	Hydromorphone-D3	20-2000	Phentermine	Phentermine-D5	100-5000
Lorazepam	Lorazepam-D4	20-1000	Pregabalin	Pregabalin-D6	20-1000
MDA	MDA-D5	100-5000	Ritalinic acid	Ritalinic acid-D10	20-1000
MDEA	MDEA-D5	100-5000	Tapentadol	Tapentadol-D3	20-1000
MDMA	MDMA-D5	100-5000	Temazepam	Temazepam-D5	20-1000
Meperidine	Meperidine-D4	20-1000	Tramadol	Tramadol-13C, D3	20-1000
Meprobamate	Meprobamate-D7	20-1000			

For research use only. Not for use in diagnostic procedures



Results


Analyte identity criteria:

- Accurate Mass of the Parent: Full-scan extracted ion chromatograms (XIC) ± 5 ppm
- Retention Time
- Isotopic Abundance Pattern
- At least two expected fragments (± 10 ppm) from the MS/MS scans
- Library Search: Custom MS/MS library built in-house using MS/MS spectra derived from calibrators



Results- MD 17500 Resolution


Internal Standard MDA-D5 is not detected in this donor sample even though it was added.





Results- MD 17500 Resolution

440130052935 #835-902 RT: 1.60-1.69 AV: 17 NL: 6.66E7
T: FTMS + p ESI Full ms [130.00-480.00]



RT 1.68 min

MDA-D5


185.1333 (± 5 ppm window)

185.1324-185.1342



Results- MD 35000 Resolution


Internal Standard MDA-D5 is not detected in this donor sample even though it was added.






Results- MD 35000 Resolution

440130052935 #649 RT: 1.69 AV: 1 NL: 1.39E7
T: FTMS + p ESI Full ms [130.00-480.00]



RT 1.68 min
MDA-D5
185.1333 (± 5 ppm window)
185.1324-185.1342

Results- MD 70000 Resolution






Results- MD 70000 Resolution

440130052935 #376-433 RT: 1.62-1.75 AV: 15 NL: 3.56E7

T: FTMS + p ESI Full ms [130.00-480.00]



RT 1.68 min

MDA-D5

185.1333 (± 5 ppm window)


185.1324-185.1342




Results- MD All resolutions

MDA-D5 185.1333 (± 5 ppm window) 185.1324-185.1342


17 500



35 000




70 000



Results- MD 17500

Internal Standard Morphine-D3 is not detected in this donor sample even though it was added.




Positive for morphine

RT 0.72 min
Morphine-D3



Results- MD 17500

150400486771 #136-280 RT: 0.68-0.86 AV: 72 NL: 1.71E6
T: FTMS + p ESI Full ms [130.00-480.00]




RT 0.75 min
Morphine-D3 (± 5 ppm window)
289.1626
289.1612-289.1640

Results- MD 35000

Positive for morphine

RT 0.72 min
Morphine-D3






Results- MD 35000

150400486771 #69-113 RT: 0.66-0.74 AV: 23 NL: 1.62E6

T: FTMS + p ESI Full ms [130.00-480.00]



RT 0.72 min

Morphine-D3 (± 5 ppm window)


289.1626

289.1612-289.1640

Results- MD 70000

Positive for morphine


RT 0.72 min
Morphine-D3





Results- MD 70000

150400486771 #50-73 RT: 0.70-0.77 AV: 12 NL: 2.06E6
T: FTMS + p ESI Full ms [130.00-480.00]




RT 0.72 min
Morphine-D3 (± 5 ppm window)
289.1626
289.1612-289.1640




Results- MD All Resolutions

Morphine-D3 (± 5 ppm window) 289.1626 (289.1612-289.1640)


17 500



35 000



70 000



Results- MD MS/MS vs ddMS²

- ▶ **Full MS OR MS/MS** : Mode where full MS spectrum scans and MS² scans are collected . Compounds included in the inclusion list have a RT time window where MS² scans are collected.

- ▶ **Full MS/ddMS²** : Mode where full MS spectrum scans and MS² scans are collected. MS/MS is performed by the software on the top ions observed in the MS spectrum. The ions are automatically selected on the basis of their signal intensity and/or scheduling time.

- ▶ DIA proved to be superior to ddMS² to collect consistently at least one MS² scan for all the analytes at the LLOQ level.

Results

Validation criteria:

- Quantitation performed using the analyte monoisotopic mass extracted from full-scan ion chromatograms (XIC)
- Specificity/Selectivity
- Linearity ($R^2 \geq 0.99$)
- Repeatability ($RSD \leq 20\%$)
- Recovery and Accuracy Standards and QCs ($RSD \leq 20\%$)
- Stability
- Matrix interferences ($\%RSD \leq 15\%$ for the internal standard-normalized matrix factor)
- Interferences from OTC drugs, blood contamination (5% v/v).

Results Validation

Precision and Accuracy				
QC	Accuracy (% Dev)		Precision (% RSD)	
	Intra (n=6)	Inter (n=18)	Intra (n=6)	Inter (n=18)
LLOQ	≤18.8	≤ 9.7	≤ 10.5	≤ 13.4
Low	≤ 17.8	≤ 6.7	≤ 5.7	≤ 13.9
Mid	≤ 6.6	≤ 3.2	≤ 7.4	≤ 7.1
High	≤ 12.9	≤ 8.3	≤ 5.4	≤ 7.6
Correlation Coefficient (r^2)	≥ 0.997			

Results Validation


Stability		
Stability	Condition	Accuracy (% Dev)
Freeze / Thaw	3 Cycles, -20°C	≤ 9.8
Room Temperature	24 hrs	≤ 15.7
Autosampler Stability	1 week, 5°C	≤ 12.0
Long – Term Storage Stability	4 week, 5°C	≤ 17.7

Results Validation

- ▶ The internal standard-normalized matrix factor for nine independent matrix sources was $\leq 14.0\%$ RSD for all the analytes.
- ▶ Assessments of specificity, carryover, impact from blood contamination (5% v/v), and OTC drugs were also performed, and found to be acceptable.

Results: Screening criteria example



Nordiazepam Retention Time and Parent Accurate Mass



Expected RT	Actual RT	RT Delta
4.83	4.82	0.02
m/z (Expected)	m/z (Apex)	m/z (Delta(ppm))
271.0633	271.0628	-1.8255


Results: Screening criteria example

Nordiazepam Isotopic Abundance Pattern






Results: Screening criteria example

Nordiazepam Isotopic Abundance Pattern




Results: Screening criteria example

Nordiazepam Fragments Accurate Mass




Results: Screening criteria example


Nordiazepam Library Match




Results



Results



Results



Conclusion

- ▶ Resolution value of 70 000 was necessary for the full scan data.
- ▶ MS/MS proved to be superior to ddMS² to collect consistently at least one MS² scan for all the analytes at the LLOQ level.
- ▶ The optimal conditions that yielded consistent hydrolysis values $\geq 80\%$ for all the commercially available glucuronide standards were incubation at 65°C for one hour, with 5000 Fishman units/mL of enzyme concentration in the master mix solution. Codeine 6-glucuronide was the analyte that took longer to hydrolyze.
- ▶ By combining the high resolution capabilities of the Thermo Scientific™ Q Exactive™ MS with the high rate of hydrolysis of the IMCSzyme it was possible to decrease sample analysis cost while increasing accuracy and the number of analytes screened and quantified on a single injection.

Acknowledgments

Thermo Fisher Scientific:

- ▶ Marketing group
- ▶ Russ Constantineau
- ▶ Susan DiPietro

