



High Resolution MS in Forensic Toxicology Screening

Osama Abu-Nimreh
CMD Sales Support Specialist
MECEC , Dubai

Screening Approaches in LC/MSMS

- Screening applications are commonly used in forensic and clinical toxicology laboratories.
 - **Targeted screening** : compound is identified and confirmed using databases and/or libraries.
 - **Unknown screening**: no databases and libraries available. Compound is identified using MS² and or MSⁿ data.
- Screening applications utilize different types of mass spectrometers
 - **Ion Traps** : MS and MSⁿ experiments. Pos/Neg switching
 - **Triple quadrupole** : 2 SRMs/analyte. Confirmation using the Ion Ratio.
 - **HRAM instruments** (OrbiTrap) : Full Scan followed by AIF for the Exactive Plus. Full Scan followed by MS² experiments for the Thermo Scientific™ Q-Exactive™ Plus. **Full Scan followed by 4 vDIA events for the Thermo Scientific™ Q-Exactive Focus.**

Screening – General Workflows

Step 1: Sample Preparation

Depending on sample type (*urine, plasma, serum, whole blood*):

- Dilution
- LLE
- SPE
- Online TurboFlow extraction
- Protein precipitation

Step 2 : Data acquisition

Different approaches:

- Ion trap
- Triple quadrupole
- Orbitrap (HRAM MS)


Step 3: Processing Data

TraceFinder

- ToxID
- HRAM screening


Mass Resolution FWHM

- Resolution**

$$R = \frac{m}{\Delta m}$$


- Quadrupole MS

$$R = \frac{400}{0.4} = 1000$$



- Orbitrap (HRAM) MS

$$R = \frac{400}{0.004} = 100000$$



How Accurate Is Your Mass?




- Mass accuracy

$$\Delta m / z = \frac{m_{meas} - m_{true}}{m_{true}} \cdot 10^6$$

- Quadrupole MS $\Delta m / z = \frac{500.1 - 500.0}{500} \cdot 10^6 = 200 ppm$


- Orbitrap MS
TOF MS $\Delta m / z = \frac{500.10314 - 500.10214}{500.10314} \cdot 10^6 = 2 ppm$

Selectivity Increases With Higher Mass Accuracy




Specificity = Resolution + Mass Accuracy


Resolution: 10k, 30k, 50k, 100k




The Industry's Leading Thermo Scientific™ MS Portfolio



Q Exactive MS - a 3D view




Q Exactive™ : Hardware




1. Ions are injected through the source
2. ...and trapped in the C-trap and squeezed into a smaller cloud
3. ...then a voltage pulse across C-trap ejects ions towards the Orbitrap
4. ...where they are trapped and detected

Orbitrap - Principle of Operation



- A short ion packet of one m/z from c- trap enters the field tangentially
- C-trap is only used as an ion storage device
- Ions are squeezed towards the central electrode by increasing voltage on the central electrode
- In the axial direction, ions are forced to move away from the narrow gap towards the wider gap near the equator. This initiates axial oscillations
- After the voltage increase stops, ion trajectories become a stable spiral

Orbitrap - Principle of Operation



- Ion m/z separation depends on
 - Frequency of harmonic oscillations and is proportional to sq root of m/z
- Three frequencies create oscillations
 - Frequency of rotation
 - Frequency of radial oscillations
 - Frequency of axial oscillations


- Resolving power is
 - Inversely proportional to the square root of m/z
 - Proportional to acquisition time

- Sensitivity is independent of acquisition speed.


■ Red rings smallest m/z ; Blue ring larger m/z ; Green ring largest m/z

Makarov A. *Anal. Chem.* 2000, 72, 1156-1162.

Principle of Orbitrap MS Operation



$$\omega_z = \sqrt{\frac{k}{m/q}}$$



Hyper-logarithmic potential distribution:
“ideal Kingdon trap”

$$U(r, z) = \frac{k}{2} \cdot \left\{ z^2 - r^2 / 2 + R_m^2 \cdot \ln(r / R_m) \right\}$$

■ Characteristic frequencies:

- Frequency of rotation ω_ϕ
- Frequency of radial oscillations ω_r
- Frequency of axial oscillations ω_z

$$\omega_\phi = \frac{\omega_z}{\sqrt{2}} \sqrt{\left(\frac{R_m}{R}\right)^2 - 1} \quad \omega_r = \omega_z \sqrt{\left(\frac{R_m}{R}\right)^2 - 2}$$

Makarov A. *Anal. Chem.* 2000, 72, 1156-1162.

Triple Quadrupole is great tool! but ?

- It is only targeted!
- Selectivity provided by tandem MS/MS (SRM transition needed)
- False positives are reality!
- Need to setup instrument (SRM) before analysis
- Realistic breakpoint is 200-300 compounds in a run
- Time consuming data processing


HRAM is a solution!

- Can perform the same level of quantitation as MS/MS
- Selectivity obtained by accurate mass measurement (only m/z needed)
- No false positives!
- No need to setup instrument (SRM) before analysis
- Unlimited number of compounds in a run – perfect for screening
- Automated data processing

Orbitrap Mass Analyzer Features

- Fundamental difference to other HRAM instruments
- Parameter measured is **frequency**, not time/voltage/current
- Resolution allows more accurate m/z determination
- Less prone to ambient conditions changes
- Usually stable within <2 ppm during several days
- No need for lock mass in “routine work”
- Small footprint
- Easy to setup

High Resolution MS Technology Race





Orbitrap technology – Workflow examples

Step 1: Sample Preparation

Step 2 : Data acquisition

Step 3: Processing Data



Q Exactive Focus - Acquisition approaches

- There are 3 approaches possible for screening:
 - **DDE** : Data Dependent Experiment. Here the system selects the more intense ions reported in the Full scan MS spectra to fragment those on MS² mode. If the ion has a low intensity it is probable that it won't be selected for MS² and therefore not confirmed by the processing software.
 - **AIF** : All Ion Fragmentation. Here the system fragments all the ions present in the MS spectra in the collision cell. Lack of specificity.
 - **DIA** : Data Independent Analysis. Here fragmentation is performed in different mass ranges. It is more specific than AIF but less specific than DDE.

Today we use the two approaches DDE and vDIA for screening purposes, we strongly suggest the vDIA approach for a better fragmentation.



Drug identification based on :
Accurate mass of the parent ion
Accurate mass of the fragment ions
Isotopic pattern
Library match
Chromatographic retention time window

3 ways of Quantitation/Screening for Routine Work

Full MS or targeted SIM/ddMS2

- Post-acquisition - extracted ion chromatograms of parent ions of interest
- Relies on high resolution for selectivity
- Useful for less complex background
- No method development/preparation needed



Experiments



Full MS/ All Ion Fragmentation – vDIA*

- Post-acquisition - extracted ion chromatograms of parent ions of interest
- Scheduled target (inclusion) list (Rt, m/z)
- Minimum method development (e.g., predefine parent ions, tr)
- Also for screening purposes



Experiments





PRM (Parallel Reaction Monitoring)

- Post-acquisition – extracted ion chromatograms of parent -> fragment transitions acquired
- Scheduled target list (Rt, m/z , collision energy)
- Most sensitive and selective even in highly complex matrices

Experiments



The Screening Method: General features



- **LC Method**

- **Mobile Phase:**

- A: 10 mM ammonium formate in 0.1% formic acid
 - B: ACN containing 0.1% formic acid

- **LC column: PFP, 150 x 2.1 mm, 5 μ m**
 - **Injection volume 20 μ L**
 - **30 mn or 15 mn Gradient**

Screening LC conditions

- **Mobile Phases**

- **A: H₂O, [NH₄]⁺ [HCOO]⁻ 2mM, 0.1% HCOOH**

(for 1L of mobile phase A use 1L of water and add 126mg of ammonium formiate and 1mL of formic acid)


- **B: [NH₄]⁺ [HCOO]⁻ 2mM, MeOH/ACN 50/50, 0.1% HCOOH, 1 % H₂O**

(for 1L of mobile phase B use 495mL of methanol, 495mL of acetonitrile, 10mL of water, and add 126mg of ammonium formiate and 1mL of formic acid)

- **UHPLC Separation**

- Accucore Phenyl Hexyl 100 x 2.1 mm, 2.6 µm

- Column Oven : 40° C



	Retention [min]	Flow [ml/min]	%B	%C	%D	Curve
1	0.000	0.500	1.0	0.0	0.0	
2	0.000	0.500	1.0	0.0	0.0	
3	1.000	0.500	1.0	0.0	0.0	
4	10.000	0.500	99.0	0.0	0.0	
5	11.500	0.500	99.0	0.0	0.0	
6	11.700	0.500	1.0	0.0	0.0	
7	17.000	0.500	1.0	0.0	0.0	

TraceFinder 4.1

- Easy to use software for all LC MS & GC MS quantitation and screening needs
 - **User security/audit trails:** Individuals or domain groups can be given different levels of access to the system and data
 - **Common confirmations in Quan and Screening workflows: Quantitate the things that you know and screen for suspects in a single method**
 - Screening to quantitation workflows for efficient method development
 - **Enhanced custom reports** with many of the same formula functions as MS Excel for calculations and conditional formatting
 - **Intelligent Sequencing** to save time and samples



**Software for
Targeted and Non-Targeted Analysis**

- Intuitive software for routine semi-quantitation and targeted screening needs in Clinical Research and Forensic Toxicology
- Customizable databases, compound confirmations, data review layouts and reporting
- Experiment specific design for SRM, Full MS – AIF, and Full MS-data dependent MS²
- Same theme as TraceFinder
- Security

Get Results

Quick, Effortless, Accurate



ToxFinder 1.0

Targeted Screening Software

© Copyright 2014 Thermo Fisher Scientific Inc.
All rights reserved. This program is protected
by copyright law and international treaties as
described in Help About.

Thermo
SCIENTIFIC

TraceFinder : Screening workflow

Creation of the master method:

- Acquisition parameters (instrumental method)
- Processing parameters (peak detection and screening criteria)
- Reporting parameters (optional)



Acquisition of the data:

- Creation of the analysis sequence



Review the data and run reports

- Data processing is automatically performed after the injections

Note: You can also process data that have already been acquired.





Targeted / Untargeted Screening workflow

DATA FILE

PROCESSING METHOD

RESULTS

Orbitrap Data




Target List


#	Library Name	Library RT	Library MS	Library MS	Library MS
1	Ch-Ethanol	1.082402	223.0734	+	104-Gal
2	Ch-Ethanol	1.082402	223.0734	+	104-Gal
3	1-Heptyl-2-ethyl-2-methyl-2-butanol	1.082402	223.0734	+	104-Gal
4	1-Heptyl-2-ethyl-2-methyl-2-butanol	1.082402	223.0734	+	104-Gal
5	1-Heptyl-2-ethyl-2-methyl-2-butanol	1.082402	223.0734	+	104-Gal
6	1-Heptyl-2-ethyl-2-methyl-2-butanol	1.082402	223.0734	+	104-Gal
7	1-Heptyl-2-ethyl-2-methyl-2-butanol	1.082402	223.0734	+	104-Gal
8	1-Heptyl-2-ethyl-2-methyl-2-butanol	1.082402	223.0734	+	104-Gal
9	1-Heptyl-2-ethyl-2-methyl-2-butanol	1.082402	223.0734	+	104-Gal
10	1-Heptyl-2-ethyl-2-methyl-2-butanol	1.082402	223.0734	+	104-Gal
11	1-Heptyl-2-ethyl-2-methyl-2-butanol	1.082402	223.0734	+	104-Gal
12	1-Heptyl-2-ethyl-2-methyl-2-butanol	1.082402	223.0734	+	104-Gal
13	1-Heptyl-2-ethyl-2-methyl-2-butanol	1.082402	223.0734	+	104-Gal
14	1-Heptyl-2-ethyl-2-methyl-2-butanol	1.082402	223.0734	+	104-Gal
15	1-Heptyl-2-ethyl-2-methyl-2-butanol	1.082402	223.0734	+	104-Gal

Searches for peaks

Compound m/z and RT




Isotope Pattern Match



Full Scan

Identification and confirmation using m/z, RT, isotope pattern

Library and Fragment Match



AIF Scan, MS² Scan

Identification and confirmation using fragment and library matching

Identification Confirmed !
Acetamidrid

Data review - flags

The table results: flag code

General flags that are the results of the individuals flags (same flags than in the Sample list)


- : Compound has been detected and confirmed -> All the individual flags are OK
- ▲ : Compound has been detected but not confirmed in the sample -> At least one of the confirmation criteria (IP or LS in this example) is in red
- : Compound has not been detected in the sample -> the identification criteria (m/z and RT in this example) are in red

Individual flags for m/z, RT, Isotopic pattern and Library search according to the criteria specified in the master method.

- : The m/z, RT, IP or LS is inside the limits specified in the master method
- : The m/z, RT, IP or LS is outside the limits specified in the master method

Flag	MZ	IP	FI	Compound Name	Formula	Confirmed	m/z (Expected)	m/z (Apex)
●	●	●	●	11-Hydroxy-delta-9-THC	C21H30O3	2 out of 2	331.22677	331.22647
●	●	●	●	6-Monoacetylmorphine	C19H21NO4	2 out of 2	328.15433	328.15414
●	●	●	●	Amisulpride	C17H27N3O4S	2 out of 2	370.1795	370.17929
●	●	●	●	Amisulpride	C20H23N	2 out of 2	278.19033	278.19037
●	●	●	●	Amisulpride	C20H25CIN2O5	2 out of 2	409.15248	409.15234
●	●	●	●	Antipyrine	C11H12N2O	2 out of 2	189.10224	189.10234
●	●	●	●	Aterolol	C14H22N2O3	2 out of 2	267.17032	267.17017
●	●	●	●	Atropine	C17H23NO3	2 out of 2	290.17507	290.17502
●	●	●	●	Benzocaine	C9H11NO2	2 out of 2	166.08625	166.08629
●	●	●	●	Benzoylcegonine	C16H19NO3	2 out of 2	290.13868	290.13867
●	●	●	●	Biperiden	C21H29NO	2 out of 2	312.23219	312.23196
●	●	●	●	Buprenorphine	C29H41NO4	2 out of 2	468.31084	468.3107


Data review – Chromatogram and spectrum




The chromatogram view

The spectrum view


Data review – isotopic pattern



Data review - Fragments



Data review – library comparison



Screening – Excel Export of a table

Thermo TraceFinder LC

Real time status | User: valerie.thibert | ? | ⚙
Waiting for Device Initialization

File View Tools Help

New
Open
Save
Print batch

Export data to
Recent Files
Exit

Target Screening

Report View

Local Method

Acquisition
Screening
Processing
Peak Detection
Reports

Acquisition
Analysis
Method Development

Data Review - 20150310

Samples
Compounds

File	MZ	IP	FI	LS	Compound Name	Lib Match Name	Formula
					Bromazepam	N/A	C14H10BrN3O
					Cetirizine	N/A	C21H25ClN2O3
					Citalopram	N/A	C20H21FN2O
					Clobazam	N/A	C16H13ClN2O2
					Clozapepam	N/A	C15H10ClN3O3
					propoxypl	N/A	C22H29NO2
					pam	N/A	C16H13ClN2O
					nhydramin	N/A	C17H21NO
					amine	N/A	C17H22N2O
					lopram	N/A	C20H21FN2O

Data Review Export

Filename: C:\TraceFinderData\32\Projects\Screening\20150310.xlsx

File Format
☒ Excel ☐ CSV

Sheet Layout
☐ Multiple Worksheets ☒ Single Worksheet

Data To Export
☒ Export selected rows (visible columns only)
☐ Export all batch data

Export Cancel

Relative Intensity

RT(min)

m/z

NL: 7.74E8
Diazepam

39

33.07111

700 800 900 1000

Report View

Reporting - 20150310

Template

Target Screening High Density Sample Report
Target Screening High Density Sample Report 2
Target Screening Summary Report
Targeted Screening _ VT

Rules

Sheet Name	Rules
Sheet1	EachSample
Sheet2	EachSample

Customize your reports

Generate a csv, excel or PDF report

Design

New

Preview

☒ PDF

☐ Excel

☐ CSV

☐ Print

Generate

Preview a report

Report Preview - Target Screening Summary Report.xlsx

Preview a report

Report Ty	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q
1	Target Screening Summary Report																
2	Batch:	20150310							Lab Name:	Default Laboratory							
3	Method:	Screening-dd_AccPheHex_20150310							Assay:	Assay name							
4	Inst. Method:	Screening-dd_AccPheHex_20150310							User:	userqex2							
5																	
6	Sample ID	Raw File			Sample Type			Vial Position			Acquisition Date						
7		Melange_std_1ugml_01			Unknown			GE2			3/10/2015 2:47:23 PM						
8																	
9	C:\Thermo\TraceFinder\3.2\General\Databases\Screening_AccPheHex_20150310.cdb																
10	Found	Confirmed	Target Name	+/-	Area	RT (Meas)	Formula	Expected m/z	Measured m/z	Delta m/z	Isotopic Pattern Score (%)	Num Isotopes Matched	Library Match	Library Score (%)	Fragments Found	Adducts	
11	0 out of 2	0 out of 2	19-Norandrosterone	+	2.58E+05	0.52	C18H28O2	277.21621	277.21664	1.57	20	1 of 2	N/A	N/A	0	Hydrogen**	
12	0 out of 2	0 out of 2	19-Noretiocholanolone	+	2.58E+05	0.52	C18H28O2	277.21621	277.21664	1.57	20	1 of 2	N/A	N/A	0	Hydrogen**	
13	0 out of 2	0 out of 2	4-Aminobiphenyl	+	1.75E+07	0.52	C12H11N	170.09643	170.09679	2.10	33	1 of 3	N/A	N/A	0	Hydrogen**	
14	0 out of 2	0 out of 2	4-Butoxyphenylacetic aci	+	3.45E+05	0.52	C12H16O3	209.11722	209.11717	-0.23	35	1 of 2	N/A	N/A	0	Hydrogen** Hydrogen-	
15	1 out of 2	1 out of 2	Alpha-Thujone	+	8.24E+06	0.52	C10H16O	153.12739	153.12767	1.83	97	2 of 2	N/A	N/A	0	Hydrogen**	
16	0 out of 2	0 out of 2	Dihydrotestosterone	+	5.97E+06	0.52	C19H30O2	291.23186	291.23187	0.04	70	2 of 3	N/A	N/A	0	Hydrogen**	

Industry Leading HRAM Library

Comprehensive HRAM Library and DB created on Thermo Scientific™ Q Exactive™ MS at R 140,000

Searchable in TraceFinder software

Consists of:

- Pesticides, Mycotoxins, Veterinary Drugs, Environmental Contaminates, PFCs
- Clin/Tox (Drugs of Abuse, Therapeutic Drugs, Poisons)
- The new spectra library will include the following: 3 ramped CE @ 20, 30, 40 eV and 2 step collision energies @ 40 with 50% and 70 with 50%
- Will contain RTs, and RRTs using the same group of ISDs for both EFS + Clin/Tox

EFS + Clin/Tox MS/MS Spectra to be available in mzCloud


Compound Groupings	Unique Entries	Total Spectra
Environmental and Food Safety	1,634	8,906
Clinical Research and Forensic Toxicology	926	4,630

Compound Classes Provided in HRAM MS/MS Libraries

COMPOUND CLASS	
Food Safety and Environmental	Forensic Toxicology
Emerging Environmental Contaminants	Drugs of Abuse
Pesticides	Natural and Industrial Toxins
Veterinary Drugs	Prescription Drugs
Mycotoxins	Performance Enhancing Drugs
Perfluorinated Compounds (PFCs)	Other Drug Monitoring Research

- A novel mass database/library of **MS/MS** and **MSⁿ** spectra (140.000 FWHM at m/z 200)
- Structural info for compounds even if they are not represented in the library through identification of **substructures**
- Multi-energy, Multi-fragment level, Multi-fragment technique
- Open consortium to establish a large public domain library which our software will link to

• <https://www.mzcloud.org/>



The screenshot shows the mzCloud website interface. At the top, there is a navigation bar with links: Home, About, Features, App, Database, Partners, and Contact. A search bar is located in the top right corner with the text "search for compounds..." and a "Search" button. Below the navigation bar, a large blue button labeled "Enter Database" is prominent. To its right, there is a section for "New mzCloud App!" with icons for Google Play and the App Store. Further right, a video player shows a screenshot of the database interface. Below these elements, there is a search section titled "Search for Compounds by Name or ID" with a text input field and an orange "Search" button. A small text block describes mzCloud as a state-of-the-art mass spectral database that assists analysts in identifying compounds in various fields like life sciences, metabolomics, pharmaceutical research, toxicology, forensic investigations, environmental analysis, food control, and various industrial applications. It mentions that mzCloud features a freely searchable collection of high resolution/accurate mass spectra using a new third generation spectra correlation algorithm. A link "read more..." is provided below this text.


6,321 (+61)
compounds

9,985 (+78)
trees

1,923,641 (+8,811)
spectra

697,276 (+1,632)
QM models

[view more
statistics](#)



PLANET ORBITRAP

Welcome Guest!

REGISTER

LOGIN

SEARCH TERMS

SEARCH

SITE

SEARCH

LIBRARY

Q

ORBI-TIPS

APPLICATIONS

PRODUCTS

LIBRARY

GRANT CENTRAL

COMMUNITY

CONTACT

What's New on PLANET ORBITRAP?

Orbitrap Fusion Lumos MS

Learn about the latest
ground-breaking innovation in
Tribrid mass spectrometry



SEE THE SCIENCE OF ASMS 2016

Translating
Innovation
into Results

Learn More



NEW LIBRARY ADDITIONS

Proteomics reveals the effects of sustained weightloss on the human plasma proteome

Philipp E Geyer, Nicolai J Wewer Albrechtsen, Stefka Tyanova, Niklas Grassl, Eva W Iepsen, Julie Lundgren, Sten Madsbad, Jens J Holst, Signe S Torekov and Matthias Mann

A Nontargeted UHPLC-HRMS Metabolomics Pipeline for Metabolite Identification: Application to Cardiac Remote Ischemic Preconditioning

Judith Kouassi Nzoughe, Cinzia Bocca, Gilles Simard, Delphine Duvion-Mabean, Juan Manuel Chao de la Plaza, Dominique

