


**Next-generation sequencing**

# Multicenter in-house evaluation of an amplicon-based NGS panel for comprehensive molecular profiling

Summary of original research article published in *Molecular Diagnosis & Therapy*, Jan 11, 2025 [1]

## Study overview


- A multicenter evaluation of the Ion Torrent™ OncoPrint™ Comprehensive Assay Plus for comprehensive genomic profiling (CGP) was performed using a cohort of pan-cancer research samples (>13 different tumor types, Figure 1).
- 193 pan-cancer research samples (125 DNA and 68 RNA samples) were analyzed to evaluate concordance of the OncoPrint Comprehensive Assay Plus with orthogonal methods.
- The research samples were pre-characterized for single nucleotide variants (SNVs), insertions and deletions (indels), copy number variations (CNVs), fusions, homologous recombination deficiency (HRD), tumor mutational burden (TMB), and microsatellite instability (MSI).
- Reproducibility of results with replicate samples was also assessed across 5 academic and clinical research centers located in Europe.
- With an easy, highly automated workflow, the OncoPrint Comprehensive Assay Plus allowed for CGP from >500 genes, including genomic signatures.



**Figure 1. Distribution of tumor types for the analytical evaluation of the OncoPrint Comprehensive Assay Plus.** CNS: central nervous system. CUP: cancer of unknown primary. GIST: gastrointestinal stromal tumor. Other: includes tumor types represented by a single sample.

## High sequencing success rate, even with challenging samples

The cohort of samples was composed of diverse tumor types with variable tumor cell content (Figure 2). Some samples contained only 10%–20% tumor cell content. Despite this heterogeneity, the OncoPrint Comprehensive Assay Plus achieved a high overall (DNA and RNA) sequencing success rate of ~94% (Table 1).



**Figure 2. Variability of tumor cell content in pan-cancer samples evaluated.**

**Table 1. Sequencing success rate of the OncoPrint Comprehensive Assay Plus.**

	Samples (n)	Success rate
DNA	140	96.6%
RNA	61	89.7%
Overall (DNA and RNA)	201	93.9%

### Accurate detection of genomic alterations

The OncoPrint Comprehensive Assay Plus demonstrated high concordance with orthogonal methods for the detection of SNVs, indels, CNVs, and fusions (Tables 2–4).

**Table 2. Concordance of results for SNV and indel detection.**

	Expected	Observed ( $\geq 5\%$ AF*)	Observed ( $< 5\%$ AF)	Concordance**
SNVs	258	241	5	95.3%
Indels	32	28	1	90.6%
<b>Total (SNVs and indels)</b>	<b>290</b>	<b>269</b>	<b>6</b>	<b>94.8%</b>

\* AF: allele fraction.

\*\* Concordance was calculated irrespective of the 5% AF limit.

**Table 3. Concordance of results for CNV detection.**

	Expected	Observed	Concordance
Copy number gain ( $\geq 6$ copy number)	44	44	100.0%
Copy number loss ( $\leq 0.5$ -fold difference)	13	11	84.6%
<b>Total CNVs</b>	<b>57</b>	<b>55</b>	<b>96.5%</b>

**Table 4. Concordance of results for fusion detection.**

		Orthogonal method	
		Positive	Negative
OncoPrint Comprehensive Assay Plus	Positive	49	1
	Negative	2	0
	<b>Concordance</b>	<b>94.2%</b>	

### Accurate detection of genomic signatures

The OncoPrint Comprehensive Assay Plus demonstrated high concordance with orthogonal methods for the detection of HRD, TMB, and MSI (Table 5).

**Table 5. Concordance of results for HRD, TMB, and MSI detection.**

	Cancer types	Samples (n)	Concordance
HRD	Ovarian and breast cancer	18	100.0%
TMB ( $\geq 10$ mutations/Mb)	Pan-cancer	32	81.3%
MSI	Pan-cancer	26	80.8%
MSI	Colorectal cancer	10	100.0%

### High reproducibility demonstrated between laboratories

Across 5 laboratories, the OncoPrint Comprehensive Assay Plus demonstrated a high level of reproducibility, with a median reproducibility for SNV and indel detection of 94.0% overall. Testing conducted on replicate samples across the laboratories also showed consistency for detecting HRD, TMB, and MSI.

## Conclusions

- Based on this study, the OncoPrint Comprehensive Assay Plus demonstrated excellent analytical performance as part of a sensitive and specific platform for CGP.
- The overall concordance between the OncoPrint Comprehensive Assay Plus and orthogonal methods was 94.8% for SNVs and indels, 96.5% for CNVs, and 94.2% for fusions.
- The OncoPrint Comprehensive Assay Plus also demonstrated high concordance with orthogonal methods for the detection of HRD, TMB, and MSI genomic signatures.
- The high concordance rates highlight the potential of the assay to advance research in the field of personalized medicine.\*

\* The OncoPrint Comprehensive Assay Plus is for research use only, and this analysis was performed as part of a retrospective clinical research study; no patient management decisions were made based on these results.

## Reference

1. Jantus-Lewintre E, Rappa A, Ruano D et al. (2025) Multicenter in-house evaluation of an amplicon-based next-generation sequencing panel for comprehensive molecular profiling. *Mol Diagn Ther* 29(2):249–261. [doi.org/10.1007/s40291-024-00766-2](https://doi.org/10.1007/s40291-024-00766-2)

 Learn more about the OncoPrint Comprehensive Assay Plus and watch the webinar at [thermofisher.com/ocaplus](https://thermofisher.com/ocaplus)