




Dioxin and PCBs with Triple Quadruple and HRMS Technology- Setting the New Standard


Elsamoul Hamdnalla
Thermo Fisher Scientific

1. What are POPs and Dioxins?
2. Why do we analyze for Dioxins?
3. How do we analyze for Dioxins?
 - Analysis
 - Screening
 - Confirmation
 - New EU legislation


What is a Dioxin?



Polychlorinated dibenzo-p-dioxin (PCDD)
"Dioxin"
75 congeners



Polychlorinated dibenzofuran (PCDF)
"Furan"
135 congeners



Polychlorinated Biphenyl (PCB)
209 congeners

**29 toxic dioxin and
PCB congeners
out of total of 419**

'Dioxins' – Complex Mixture Of Compounds

- Polychlorinated dibenzodioxins –PCDDs (75)
 - Toxic congeners: 7
- Polychlorinated dibenzofurans - PCDFs (135)
 - Toxic Congeners: 10
- Polychlorinated biphenyls
 - Toxic congeners: 12

Compound Group	Mono	Di	Tri	Tetra	Penta	Hexa	Hepta	Octa	Nona	Deca	Total	Toxic
PCB	3	12	24	42	46	42	24	12	3	1	209	12
PCDD	2	10	14	22	14	10	2	1			75	7
PCDF	4	16	28	38	28	16	4	1			135	10
PBB	3	12	24	42	46	42	24	12	3	1	209	
PBDE	3	12	24	42	46	42	24	12	3	1	209	
PBDD	2	10	14	22	14	10	2	1			75	?
PBDF	4	16	28	38	28	16	4	1			135	?
PBCDD	0	14	84	254	420	452	252	74			1550	?
PBCDF	0	28	168	496	840	880	504	134			3050	?




- Viktor Yushchenko, President of Ukraine (2005-2010)
- Poisoned in 2004

Why Analyze for Dioxins?

UNEP definition of Persistent Organic Pollutants (POPs)

“...chemical substances that **PERSIST** in the environment, **BIOACCUMULATE** through the food web, and pose a **RISK** of causing adverse effects to **HUMAN HEALTH** and the environment. With the evidence of long-range transport of these substances to regions where they have never been used or produced and the consequent threats they pose to the environment of the whole globe, the international community has now, at several occasions called for urgent global actions to **REDUCE AND ELIMINATE RELEASES** of these chemicals.”



Dioxins at Trace Levels are Highly Toxic

- Tolerable Daily Intake 1 pg/kg/d TEQ
- Amounts typically detected:
 - Food PPT
 - Water PPQ
 - Human Serum PPQ
 - Air fg/m³

Toxicity of Dioxins


- Acute toxicity
 - Chloracne
 - Carcinogenic
 - Estrogenic
- President Yushchenko (Ukraine)
 - Suffered from Dioxin intoxication during his election in 2004
 - TCDD 10-12 $\mu\text{g/L}$ in blood
 - 100 **ppb** - 50,000x higher than normal concentrations



- “Dioxin so pure it was made in a laboratory”

Basic Concepts of Toxic Equivalency Factor (TEF)

- Dioxin-like toxicity is **cumulative**
 - So we can compare toxicology of various dioxin-like compounds and simplify risk assessment and regulatory control.
- All other compounds with dioxin-like toxicity are compared to 2378 TCDD and given a relative TEF



- Toxic Equivalency Factor (TEF) of 2378-TCDD = 1

Ever Changing TEF Concept

- TEF factors are continuously revised
- TEFs are species dependent
 - Human
 - Fish
 - Avian
 - Mammalian
- Generally we are interested in human exposure
- EU guidelines for food analysis use WHO2005 TEF factors

Congener	I-TEF	WHO1998-TEF	WHO2005-TEF
<i>Polychlorinated dibenzo-p-dioxins</i>			
2378-CI4DD	1	1	1
12378-CI5DD	0.5	1	1
123478-CI6DD	0.1	0.1	0.1
123678-CI6DD	0.1	0.1	0.1
123789-CI6DD	0.1	0.1	0.1
1234678-CI7DD	0.01	0.01	0.01
CI8DD	0.001	0.0001	0.0003
<i>Polychlorinated dibenzofurans</i>			
2378-CI4DF	0.1	0.1	0.1
12378-CI5DF	0.05	0.05	0.03
23478-CI5DF	0.5	0.5	0.3
123478-CI6DF	0.1	0.1	0.1
123678-CI6DF	0.1	0.1	0.1
123789-CI6DF	0.1	0.1	0.1
234678-CI6DF	0.1	0.1	0.1
1234678-CI7DF	0.01	0.01	0.01
1234789-CI7DF	0.01	0.01	0.01
CI8DF	0.001	0.0001	0.0003
<i>Non-ortho polychlorinated biphenyls</i>			
PCB-77	-	0.0001	0.0001
PCB-81	-	0.0001	0.0003
PCB-126	-	0.1	0.1
PCB-169	-	0.01	0.03
<i>Mono-ortho polychlorinated biphenyls</i>			
PCB-105	-	0.0001	0.00003
PCB 114	-	0.0005	0.00003
PCB-118	-	0.0001	0.00003
PCB-123	-	0.0001	0.00003
PCB-156	-	0.0005	0.00003
PCB-157	-	0.0005	0.00003
PCB-167	-	0.00001	0.00003
PCB-189	-	0.0001	0.00003

“Dioxin Incidents” – Economical Impact

Almost all food contamination comes from contaminated animal feed

- 1998 EU Milk from citrus pulp pellet feed
- **1999 Belgium Eggs and poultry – Total cost estimated € 625M**
- 1999 EU Clay and zeolithes for feed
- 2000 EU Choline chloride
- 2002 UK and Eire Feed supplement “Carbosan Copper”
- 2004 EU Potato pulp
- 2005 EU Hydrochloric acid
- 2007 India Guar Gum thickener
- **2008 Ireland Pork – Total cost estimated € 150M**
- 2008 Italy Mozzarella Cheese
- **2010 Germany Eggs & pork – 4700 farms closed.
Current cost estimated > € 300M**
- 2012 Germany - More poultry farms closed (local contamination)
- **2013 ????????**




- All exports to/from EU are subject to the same legislation

Who Tests For Dioxins In Food?

1. Government regulatory laboratories working in food safety and environmental
 - e.g. CVUA Frieburg, CVUA Muenster, FERA, CDC, EPA, FDA
2. Contract Testing Organisations - food safety and environmental
 - e.g. Eurofins, SGS, ALS
3. Food conglomerates working for food safety and brand protection
 - e.g. Kraft, Cargil, Coca-Cola, Thames Water

Bloomberg Businessweek Friday April 15, 2011

Get our new FREE iPad app now | Home | Finance | Technology | Innovation



Bloomberg

Pigs Don't Fly to China as Germany Struggles to Contain Crisis


January 12, 2011, 6:39 AM EST
By Patrick Donahue

Jan. 12 (Bloomberg) -- German authorities contain a widening tainted-food scandal temporarily halted imports of German pork products in response to the discovery of dioxin.

China's General Administration of Quality Inspection and Quarantine announced it its website. The move followed an order northern German state of Lower Saxony of 140 pigs at a farm after tests found the levels of dioxin, a possible carcinogen, e allowable limits.

Contaminated pork may have made its v supply, since pigs from the local farm in Verden may have been exposed to dioxin in livestock feed, a district spokesman sai the European Union's largest pig breede EU statistics agency Eurostat.

Italy defends mozzarella, Japan blocks imports



Recommend

Share on LinkedIn

Share on Facebook

Digg

Related Topics
Health »

PICTURES

SLIDESHOW

Best of the w
Our top pictures fro past week. Slides

China's richest vi

Entertainers of th

Environment wat

The Korean War

Pictures of the ye Sports


Rank	Previous Rank	Brand	Country of Origin	Sector	Brand Value (\$bn)
1	1	Coca-Cola	United States	Beverages	70,452
6	6	McDonald's	United States	Restaurants	33,578
17	20	Apple	United States	Electronics	21,143
35	34	Kellogg's	United States	FMCG	11,041
46	48	Heinz	United States	FMCG	7,534
23	23	Pepsi	United States	Beverages	14,961
57	58	Nestle	Switzerland	FMCG	6,548
81	92	Shell	Netherlands	Energy	4,003
90	NEW	3M	United States	Diversified	3,586
99	100	Campbell's	United States	FMCG	3,241

How To Analyze Dioxins - Complete POPs Solutions

Step	Requirement	Solution
Laboratory PPE	Fume hoods, lab coats, gloves, eye protection etc.	Fisher Scientific
Laboratory Reagents	Alumina, Silica, Solvents, Acids	Fisher Scientific
Laboratory Apparatus	Freeze driers, ovens, balances, pipettes, vials, columns	Fisher Scientific
Sample Preparation	Extraction	Thermo Scientific™ Dionex™ ASE350 Accelerated Solvent Extraction
	Extract Fractionation / Purification	CAPE Technologies Sample Preparation Kits*
	Extract Concentration	Thermo Scientific™ Dionex™ Rocket™ Evaporator
Sample Analysis	Chromatography	Thermo Scientific™ Trace 1310 GC, TriPlus RSH & Trace Gold GC/MS columns and consumables
Sample Analysis	GC-HRMS or GC-MS/MS	Thermo Scientific™ DFS™ GC-HRMS or Thermo Scientific™ TSQ™ 8000 Evo GC-MS/MS
Data Handling & Reporting	Dioxin Data processing, QA/QC checking & reporting	Thermo Scientific™ TargetQua™ n 3 Thermo Scientific™ Xcalibur™
	LIMS	Thermo Scientific™ Nautilus™ LIMS

- Thermo Fisher Scientific supports the entire dioxin analysis workflow

Typical Dioxin Workflow (EPA1613)



Instrumental Analysis: Instrument Characteristics Required For Pops Analysis

1. Chromatography

- Must be able to separate POP compounds & congeners

2. Sensitivity

- Compounds typically ppb to ppt concentration in environmental samples

3. Selectivity

- Must be halogen selective or have resolve analyte from matrix

4. Low price

- Environmental market is very cost sensitive

5. Matrix tolerance

- Must be capable of prolonged analysis with

6. Conformity to standard methods for POPs

- EPA, EN, MCERTS, JIS etc...

- HRGC

- GC-ECD, GC-MS, GC-MS/MS and GC-HRMS

- GC-MS/MS and GC-HRMS

- GC-ECD, GC-MS and GC-MS/MS

- GC-MS/MS and GC-HRMS

- **GC-HRMS and GC-MS/MS*** (*EU)
(for defensible confirmation)

Analytical Selectivity



TSQ 8000 Evo

DFS

GC-MS/MS

GC-HRMS

Triple Quadrupole MS

Magnetic Sector HRMS

High performance, easy to use MS/MS for non-experts

High-Resolution Full Scan and SIM

Target Analysis in Complex Matrix

Dioxins, PCBs, Other POPs

Pesticides in Food, Drugs in body fluid

EPA 1613, 8270A, 8290

ThermoFisher SCIENTIFIC Search All [dropdown] [search icon] [Contact Us] | Sign In | Quick Order [shopping cart icon]

Recently Viewed Items [dropdown]

Applications & Techniques | Shop All Products | Services & Support | About Us | Cloud


Home > Industrial & Applied Science > Environmental > Water Analysis > Wastewater Analysis > Instruments for Persistent Organic Pollutant (POPs) Analysis

Instruments for Persistent Organic Pollutant (POPs) Analysis

Wastewater Analysis

Instruments for Persistent Organic Pollutant (POPs) Analysis

Automated Discrete Analyzers
Colorimeters




Thermo Fisher Scientific has a strong commitment to supporting laboratories charged with the task of researching and monitoring POPs. Our commitment ensures high productivity, added value solutions for some of the most challenging POPs determinations.

Resources

- Power Plant Water Analysis Products
- Drinking Water Quality Analysis Process Products
- ROSS pH Electrodes Brochure
- SampleManager LIMS™ product brochure
- POPs Center of Excellence

Support

Unity Lab Services—Access to our expert service team, information on support plans, parts and consumables resources, and other instrument services.



- Two sequences were analysed on the system
 1. PCDD/Fs consisting of:-
 - CSL, CS1→CS6 EN:1948 standards in duplicate
 - Samples (fish and egg) bracketed by blank and CSL injections
 - Total sequence length 103 injection over ~4 days
 2. dl-PCBs consisting of: -
 - CS1→CS6 WP-CVS standards in duplicate
 - Samples (fish and egg) bracketed by CS4/100 diluted standards (up to 42 injections)
 - Serial dilutions of WP-CS4 standard (/100, /200, /400) to establish LOQ
 - Total sequence length 72 injections

GC and MS Conditions

TRACE 1310 GC Parameters

Injection Volume (µL):	3
Liner:	LinerGOLD™ double taper (P/N: 453A1345-UI)
Inlet (°C):	280
Carrier Gas, (mL/min):	He, 1.2
Inlet Mode,	Splitless (split flow 120mL/min after 2 min)
Oven Temperature Program:	
Temperature 1 (°C):	170
Hold Time (min):	0
Temperature 2 (°C):	250
Rate (°C/min):	15.4
Hold Time (min):	0
Temperature 3 (°C):	285
Rate (°C/min):	2.5
Hold Time (min):	0
Temperature 3 (°C):	320
Rate (°C/min):	10
Hold Time (min):	15
Total Run Time (min):	39.7

TSQ 8000 Evo Mass Spectrometer Parameters

Transfer Line (°C):	300
Ionization Type:	EI
Ion Source(°C):	350
Electron Energy (eV):	40
Acquisition Mode:	Timed SRM with Dwell Time Prioritization
Acquired Masses (Da):	Optimised
Collision Energy (V):	Optimised

Calibration Data PCDD/Fs

EN:1948 CSL, CS1-CS6 PCDD/Fs

Congener (native)	R ²	Mean RF	Std.Dev	RSD(%)
2378-TCDF	0.9999	0.999	0.029	2.9%
2378-TCDD	0.9998	1.131	0.048	4.3%
12378-PeCDF	0.9996	1.046	0.020	1.9%
23478-PeCDF	0.9999	1.086	0.026	2.4%
12378-PeCDD	0.9996	1.121	0.043	3.9%
123478-HxCDF	0.9998	1.153	0.025	2.1%
123678-HxCDF	1.0000	1.152	0.034	3.0%
234678-HxCDF	0.9997	1.117	0.025	2.2%
123478-HxCDD	0.9996	1.175	0.057	4.9%
123678-HxCDD	0.9995	1.087	0.029	2.7%
123789-HxCDD	0.9986	1.062	0.039	3.7%
123789-HxCDF	0.9990	1.058	0.035	3.3%
1234678-HpCDF	0.9999	1.078	0.022	2.1%
1234678-HpCDD	0.9998	1.070	0.034	3.2%
1234789-HpCDF	0.9999	1.112	0.023	2.1%
OCDD	0.9999	1.082	0.035	3.2%
OCDF	0.9998	1.518	0.038	2.5%


	1948CSL (pg/µl)	1948CS1 (pg/µl)	1948CS2 (pg/µl)	1948CS3 (pg/µl)	1948CS4 (pg/µl)	1948CS5 (pg/µl)	1948CS6 (pg/µl)
NATIVE PCDDs & PCDFs							
2,3,7,8-Tetrachlorodibenzo-p-dioxin	0.04	0.2	0.8	4	16	80	320
1,2,3,7,8-Pentachlorodibenzo-p-dioxin	0.08	0.4	1.6	8	32	160	640
1,2,3,4,7,8-Hexachlorodibenzo-p-dioxin	0.08	0.4	1.6	8	32	160	640
1,2,3,6,7,8-Hexachlorodibenzo-p-dioxin	0.08	0.4	1.6	8	32	160	640
1,2,3,7,8,9-Hexachlorodibenzo-p-dioxin	0.08	0.4	1.6	8	32	160	640
1,2,3,4,6,7,8-Heptachlorodibenzo-p-dioxin	0.16	0.8	3.2	16	64	320	1280
Octachlorodibenzo-p-dioxin	0.16	0.8	3.2	16	64	320	1280
NATIVE PCDFs							
2,3,7,8-Tetrachlorodibenzofuran	0.04	0.2	0.8	4	16	80	320
1,2,3,7,8-Pentachlorodibenzofuran	0.08	0.4	1.6	8	32	160	640
2,3,4,7,8-Pentachlorodibenzofuran	0.08	0.4	1.6	8	32	160	640
1,2,3,4,7,8-Hexachlorodibenzofuran	0.08	0.4	1.6	8	32	160	640
1,2,3,6,7,8-Hexachlorodibenzofuran	0.08	0.4	1.6	8	32	160	640
1,2,3,7,8,9-Hexachlorodibenzofuran	0.08	0.4	1.6	8	32	160	640
2,3,4,6,7,8-Hexachlorodibenzofuran	0.08	0.4	1.6	8	32	160	640
1,2,3,4,6,7,8-Heptachlorodibenzofuran	0.16	0.8	3.2	16	64	320	1280
1,2,3,4,7,8,9-Heptachlorodibenzofuran	0.16	0.8	3.2	16	64	320	1280
Octachlorodibenzofuran	0.16	0.8	3.2	16	64	320	1280
SAMPLING STANDARDS							
1,2,3,7,8-Pentachloro[¹⁴ C] ₅ dibenzofuran	16	16	16	16	16	16	16
1,2,3,7,8,9-Hexachloro[¹⁴ C] ₆ dibenzofuran	16	16	16	16	16	16	16
1,2,3,4,7,8,9-Heptachloro[¹⁴ C] ₇ dibenzofuran	32	32	32	32	32	32	32
EXTRACTION STANDARDS							
2,3,7,8-Tetrachloro[¹⁴ C] ₄ dibenzo-p-dioxin	16	16	16	16	16	16	16
1,2,3,7,8-Pentachloro[¹⁴ C] ₅ dibenzo-p-dioxin	16	16	16	16	16	16	16
1,2,3,4,7,8-Hexachloro[¹⁴ C] ₆ dibenzo-p-dioxin	16	16	16	16	16	16	16
1,2,3,6,7,8-Hexachloro[¹⁴ C] ₆ dibenzo-p-dioxin	16	16	16	16	16	16	16
1,2,3,4,6,7,8-Heptachloro[¹⁴ C] ₇ dibenzo-p-dioxin	32	32	32	32	32	32	32
Octachloro[¹⁴ C] ₈ dibenzo-p-dioxin	32	32	32	32	32	32	32
EXTRACTION STANDARDS							
2,3,7,8-Tetrachloro[¹⁴ C] ₄ dibenzofuran	16	16	16	16	16	16	16
2,3,4,7,8-Pentachloro[¹⁴ C] ₅ dibenzofuran	16	16	16	16	16	16	16
1,2,3,4,7,8-Hexachloro[¹⁴ C] ₆ dibenzofuran	16	16	16	16	16	16	16
1,2,3,6,7,8-Hexachloro[¹⁴ C] ₆ dibenzofuran	16	16	16	16	16	16	16
2,3,4,6,7,8-Hexachloro[¹⁴ C] ₆ dibenzofuran	16	16	16	16	16	16	16
1,2,3,4,6,7,8-Heptachloro[¹⁴ C] ₇ dibenzofuran	32	32	32	32	32	32	32
Octachloro[¹⁴ C] ₈ dibenzofuran	32	32	32	32	32	32	32
SYRINGE STANDARDS							
1,2,3,4-Tetrachloro[¹⁴ C] ₄ dibenzo-p-dioxin	16	16	16	16	16	16	16
1,2,3,7,8,9-Hexachloro[¹⁴ C] ₆ dibenzo-p-dioxin	16	16	16	16	16	16	16

WP-CVS (WP-CS1 to WP-CS6) dl-PCBs

Congener (native)	R ²	Mean RF	Std.Dev	RSD(%)
PCB -81- tetrachlorobiphenyl	0.9998	0.998	0.025	2.5%
PCB -77-tetrachlorobiphenyl	0.9999	1.032	0.027	2.6%
PCB -123- pentachlorobiphenyl	1.0000	0.937	0.011	1.2%
PCB -118 pentachlorobiphenyl	1.0000	0.993	0.016	1.6%
PCB -114 pentachlorobiphenyl	1.0000	1.037	0.027	2.6%
PCB -105 pentachlorobiphenyl	1.0000	0.956	0.020	2.1%
PCB - 126-pentachlorobiphenyl	0.9997	0.999	0.049	4.9%
PCB -167- hexachlorobiphenyl	0.9999	1.052	0.020	1.9%
PCB -156- hexachlorobiphenyl	0.9999	1.064	0.027	2.6%
PCB -157- hexachlorobiphenyl	0.9998	1.020	0.022	2.2%
PCB -169- hexachlorobiphenyl	0.9998	1.014	0.025	2.4%
PCB -189- heptachlorobiphenyl	1.0000	1.088	0.019	1.7%

NATIVE PCBs	IUPAC	WP-CS1 (ng/ml)	WP-CS2 (ng/ml)	WP-CS3 (ng/ml)	WP-CS4 (ng/ml)	WP-CS5 (ng/ml)	WP-CS6 (ng/ml)	WP-CS7 (ng/ml)
NATIVE PCBs								
3,3',4,4'-Tetrachlorobiphenyl	77	0.1	0.5	2.0	10	40	200	800
3,4,4',5-Tetrachlorobiphenyl	81	0.1	0.5	2.0	10	40	200	800
2,3,3',4,4'-Pentachlorobiphenyl	105	0.1	0.5	2.0	10	40	200	800
2,3,4,4',5-Pentachlorobiphenyl	114	0.1	0.5	2.0	10	40	200	800
2,3',4,4',5-Pentachlorobiphenyl	118	0.1	0.5	2.0	10	40	200	800
2',3,4,4',5-Pentachlorobiphenyl	123	0.1	0.5	2.0	10	40	200	800
3,3',4,4',5-Pentachlorobiphenyl	126	0.1	0.5	2.0	10	40	200	800
2,3,3',4,4',5-Hexachlorobiphenyl	156	0.1	0.5	2.0	10	40	200	800
2,3,3',4,4',5'-Hexachlorobiphenyl	157	0.1	0.5	2.0	10	40	200	800
2,3',4,4',5,5'-Hexachlorobiphenyl	167	0.1	0.5	2.0	10	40	200	800
3,3',4,4',5,5'-Hexachlorobiphenyl	169	0.1	0.5	2.0	10	40	200	800
2,3,3',4,4',5,5'-Heptachlorobiphenyl	189	0.1	0.5	2.0	10	40	200	800
MASS-LABELLED PCBs								
3,3',4,4'-Tetrachloro[¹³ C ₁₂]biphenyl	77L	50	50	50	50	50	50	50
3,4,4',5-Tetrachloro[¹³ C ₁₁]biphenyl	81L	50	50	50	50	50	50	50
2,3,3',4,4'-Pentachloro[¹³ C ₁₁]biphenyl	105L	50	50	50	50	50	50	50
2,3,4,4',5-Pentachloro[¹³ C ₁₁]biphenyl	114L	50	50	50	50	50	50	50
2,3',4,4',5-Pentachloro[¹³ C ₁₁]biphenyl	118L	50	50	50	50	50	50	50
2',3,4,4',5-Pentachloro[¹³ C ₁₁]biphenyl	123L	50	50	50	50	50	50	50
3,3',4,4',5-Pentachloro[¹³ C ₁₁]biphenyl	126L	50	50	50	50	50	50	50
2,3,3',4,4',5-Hexachloro[¹³ C ₁₁]biphenyl	156L	50	50	50	50	50	50	50
2,3,3',4,4',5'-Hexachloro[¹³ C ₁₁]biphenyl	157L	50	50	50	50	50	50	50
2,3',4,4',5,5'-Hexachloro[¹³ C ₁₁]biphenyl	167L	50	50	50	50	50	50	50
3,3',4,4',5,5'-Hexachloro[¹³ C ₁₁]biphenyl	169L	50	50	50	50	50	50	50
2,3,3',4,4',5,5'-Heptachloro[¹³ C ₁₁]biphenyl	189L	50	50	50	50	50	50	50
INTERNAL STANDARDS: MASS-LABELLED PCBs								
2,3',4',5-Tetrachloro[¹³ C ₁₁]biphenyl	70L	50	50	50	50	50	50	50
2,3,3',5,5'-Pentachloro[¹³ C ₁₁]biphenyl	111L	50	50	50	50	50	50	50
2,2',3,4,4',5'-Hexachloro[¹³ C ₁₁]biphenyl	138L	50	50	50	50	50	50	50
2,2',3,3',4,4',5-Heptachloro[¹³ C ₁₁]biphenyl	170L	50	50	50	50	50	50	50

Ion Ratio Stability



Ion ratio stability for all 17 PCDD/F congeners over 103 consecutive sample/standard injections [n=16 EN:1948 CSL (40fg/ μ L 2,3,7,8-TCDD)]

Ion Ratio Stability

Congener	Ion ratio		Deviation	StDev	RSD(%)
	Theoretical	Calculated Mean			
2378-TCDF	0.96	0.93	-3%	0.063	6.7%
2378-TCDD	0.96	0.95	-2%	0.081	8.5%
12378-PeCDF	0.78	0.82	5%	0.047	5.8%
23478-PeCDF	0.78	0.81	3%	0.054	6.6%
12378-PeCDD	0.78	0.80	2%	0.067	8.4%
123478-HxCDF	0.63	0.62	-2%	0.039	6.3%
123678-HxCDF	0.63	0.63	0%	0.052	8.2%
234678-HxCDF	0.63	0.65	4%	0.039	6.0%
123478-HxCDD	0.64	0.63	-1%	0.054	8.6%
123678-HxCDD	0.63	0.64	1%	0.045	7.0%
123789-HxCDD	0.63	0.63	0%	0.047	7.5%
123789-HxCDF	0.63	0.65	3%	0.041	6.3%
1234678-HpCDF	0.80	0.78	-3%	0.056	7.3%
1234678-HpCDD	0.80	0.79	-1%	0.033	4.2%
1234789-HpCDF	0.80	0.77	-3%	0.034	4.3%
OCDD	0.96	0.94	-2%	0.079	8.4%
OCDF	0.96	0.97	1%	0.068	7.0%

Ion ratio stability over 103 consecutive sample/standard injections (n=16), CSL ran after every 4 matrix sample injections.

Calculated Amount

Congener	Concentration (pg/ μ L)		Deviation	StDev	RSD(%)
	Nominal	Calculated Mean			
2378-TCDF	0.040	0.043	9%	0.003	7.1%
2378-TCDD	0.040	0.039	-2%	0.003	7.2%
12378-PeCDF	0.080	0.080	0%	0.003	3.8%
23478-PeCDF	0.080	0.079	-1%	0.005	5.9%
12378-PeCDD	0.080	0.077	-4%	0.006	7.8%
123478-HxCDF	0.080	0.084	5%	0.004	5.1%
123678-HxCDF	0.080	0.083	4%	0.004	4.7%
234678-HxCDF	0.080	0.079	-1%	0.005	5.8%
123478-HxCDD	0.080	0.078	-2%	0.007	8.4%
123678-HxCDD	0.080	0.083	4%	0.004	5.0%
123789-HxCDD	0.080	0.082	2%	0.007	8.0%
123789-HxCDF	0.080	0.081	2%	0.008	9.2%
1234678-HpCDF	0.160	0.171	7%	0.009	5.1%
1234678-HpCDD	0.160	0.178	11%	0.007	4.1%
1234789-HpCDF	0.160	0.164	3%	0.010	6.4%
OCDD	0.160	0.176	10%	0.014	7.8%
OCDF	0.160	0.176	10%	0.010	5.5%

Calculated concentration over 103 consecutive sample/standard injections (n=16), CSL ran after every 4 matrix sample injections.

Limit of Confirmation (LOQc)





Figure 2. Practically determined LOQc for PCDD/Fs

Serially diluted CS3 standards from CS3/100-CS3/1000

Ion ratio stability within $\pm 15\%$ tolerance

	1948CS3 (pg/ μ l)
NATIVE PCDDs & PCDFs	
2,3,7,8-Tetrachlorodibenzo-p-dioxin	4
1,2,3,7,8-Pentachlorodibenzo-p-dioxin	8
1,2,3,4,7,8-Hexachlorodibenzo-p-dioxin	8
1,2,3,6,7,8-Hexachlorodibenzo-p-dioxin	8
1,2,3,7,8,9-Hexachlorodibenzo-p-dioxin	8
1,2,3,4,6,7,8-Heptachlorodibenzo-p-dioxin	16
Octachlorodibenzo-p-dioxin	16
2,3,7,8-Tetrachlorodibenzofuran	4
1,2,3,7,8-Pentachlorodibenzofuran	8
2,3,4,7,8-Pentachlorodibenzofuran	8
1,2,3,4,7,8-Hexachlorodibenzofuran	8
1,2,3,6,7,8-Hexachlorodibenzofuran	8
1,2,3,7,8,9-Hexachlorodibenzofuran	8
2,3,4,6,7,8-Hexachlorodibenzofuran	8
1,2,3,4,6,7,8-Heptachlorodibenzofuran	16
1,2,3,4,7,8,9-Heptachlorodibenzofuran	16
Octachlorodibenzofuran	16
SAMPLING STANDARDS	
1,2,3,7,8-Pentachloro[$^{13}\text{C}_{12}$]dibenzofuran	16
1,2,3,7,8,9-Hexachloro[$^{13}\text{C}_{12}$]dibenzofuran	16
1,2,3,4,7,8,9-Heptachloro[$^{13}\text{C}_{12}$]dibenzofuran	32
EXTRACTION STANDARDS	
2,3,7,8-Tetrachloro[$^{13}\text{C}_{12}$]dibenzo-p-dioxin	16
1,2,3,7,8-Pentachloro[$^{13}\text{C}_{12}$]dibenzo-p-dioxin	16
1,2,3,4,7,8-Hexachloro[$^{13}\text{C}_{12}$]dibenzo-p-dioxin	16
1,2,3,6,7,8-Hexachloro[$^{13}\text{C}_{12}$]dibenzo-p-dioxin	16
1,2,3,4,6,7,8-Heptachloro[$^{13}\text{C}_{12}$]dibenzo-p-dioxin	32
Octachloro[$^{13}\text{C}_{12}$]dibenzo-p-dioxin	32
2,3,7,8-Tetrachloro[$^{13}\text{C}_{12}$]dibenzofuran	16
2,3,4,7,8-Pentachloro[$^{13}\text{C}_{12}$]dibenzofuran	16
1,2,3,4,7,8-Hexachloro[$^{13}\text{C}_{12}$]dibenzofuran	16
1,2,3,6,7,8-Hexachloro[$^{13}\text{C}_{12}$]dibenzofuran	16
2,3,4,6,7,8-Hexachloro[$^{13}\text{C}_{12}$]dibenzofuran	16
1,2,3,4,6,7,8-Heptachloro[$^{13}\text{C}_{12}$]dibenzofuran	32
Octachloro[$^{13}\text{C}_{12}$]dibenzofuran	32
SYRINGE STANDARDS	
1,2,3,4-Tetrachloro[$^{13}\text{C}_{12}$]dibenzo-p-dioxin	16
1,2,3,7,8,9-Hexachloro[$^{13}\text{C}_{12}$]dibenzo-p-dioxin	16




Chromatography – Lowest Calibration Standard (EN:1948 CSL)



1948CSL (pg/μl)	
NATIVE PCDDs & PCDFs	
2,3,7,8-Tetrachlorodibenzo-p-dioxin	0.04
1,2,3,7,8-Pentachlorodibenzo-p-dioxin	0.08
1,2,3,4,7,8-Hexachlorodibenzo-p-dioxin	0.08
1,2,3,6,7,8-Hexachlorodibenzo-p-dioxin	0.08
1,2,3,7,8,9-Hexachlorodibenzo-p-dioxin	0.08
1,2,3,4,6,7,8-Heptachlorodibenzo-p-dioxin	0.16
Octachlorodibenzo-p-dioxin	0.16
EXTRACTION STANDARDS	
2,3,7,8-Tetrachlorodibenzofuran	0.04
1,2,3,7,8-Pentachlorodibenzofuran	0.08
2,3,4,7,8-Pentachlorodibenzofuran	0.08
1,2,3,4,7,8-Hexachlorodibenzofuran	0.08
1,2,3,6,7,8-Hexachlorodibenzofuran	0.08
1,2,3,7,8,9-Hexachlorodibenzofuran	0.08
2,3,4,6,7,8-Hexachlorodibenzofuran	0.08
1,2,3,4,6,7,8-Heptachlorodibenzofuran	0.16
1,2,3,4,7,8,9-Heptachlorodibenzofuran	0.16
Octachlorodibenzofuran	0.16
SAMPLING STANDARDS	
1,2,3,7,8-Pentachloro[¹³ C ₁₂]dibenzofuran	16
1,2,3,7,8,9-Hexachloro[¹³ C ₁₂]dibenzofuran	16
1,2,3,4,7,8,9-Heptachloro[¹³ C ₁₂]dibenzofuran	32
EXTRACTION STANDARDS	
2,3,7,8-Tetrachloro[¹³ C ₁₂]dibenzo-p-dioxin	16
1,2,3,7,8-Pentachloro[¹³ C ₁₂]dibenzo-p-dioxin	16
1,2,3,4,7,8-Hexachloro[¹³ C ₁₂]dibenzo-p-dioxin	16
1,2,3,6,7,8-Hexachloro[¹³ C ₁₂]dibenzo-p-dioxin	16
1,2,3,4,6,7,8-Heptachloro[¹³ C ₁₂]dibenzo-p-dioxin	32
Octachloro[¹³ C ₁₂]dibenzo-p-dioxin	32
EXTRACTION STANDARDS	
2,3,7,8-Tetrachloro[¹³ C ₁₂]dibenzofuran	16
2,3,4,7,8-Pentachloro[¹³ C ₁₂]dibenzofuran	16
1,2,3,4,7,8-Hexachloro[¹³ C ₁₂]dibenzofuran	16
1,2,3,6,7,8-Hexachloro[¹³ C ₁₂]dibenzofuran	16
2,3,4,6,7,8-Hexachloro[¹³ C ₁₂]dibenzofuran	16
1,2,3,4,6,7,8-Heptachloro[¹³ C ₁₂]dibenzofuran	32
Octachloro[¹³ C ₁₂]dibenzofuran	32
SYRINGE STANDARDS	
1,2,3,4-Tetrachloro[¹³ C ₁₂]dibenzo-p-dioxin	16
1,2,3,7,8,9-Hexachloro[¹³ C ₁₂]dibenzo-p-dioxin	16

Natives only

Sample – Selected Chromatography And Calculations



123789-HxCDF – “below” LOQc


12378-PeCDD – 0.162 pg/g WHO-TEQ

2378-TCDD – 0.193 pg/g WHO-TEQ

Peak Name	Retention Time min	Amount pg/μL	LOQc pg/μL	WHO-TEF (2005)	WHO-TEQ(2005) pg/g	Peak Confirmation
2378-TCDF	18.642	12.8690	0.02	0.1	0.512	Valid
2378-TCDD	19.103	0.4864	0.02	1	0.193	Valid
12378-PeCDF	21.977	0.7828	0.048	0.03	0.009	Valid
23478-PeCDF	22.956	2.1316	0.04	0.3	0.254	Valid
12378-PeCDD	23.117	0.4070	0.072	1	0.162	Valid
123478-HxCDF	25.356	0.1276	0.04	0.1	0.005	Valid
123678-HxCDF	25.456	0.1242	0.048	0.1	0.005	Valid
234678-HxCDF	26.056	0.1160	0.048	0.1	0.005	Valid
123478-HxCDD	26.118	0.0938	0.064	0.1	0.004	Valid
123678-HxCDD	26.207	0.2346	0.064	0.1	0.009	Valid
123789-HxCDD	26.507	0.0825	0.072	0.1	0.003	Valid
123789-HxCDF	26.878	0.0301	0.064	0.1	0.001	nd<LOQc
1234678-HpCDF	28.396	0.1036	0.112	0.01	0.000	nd<LOQc
1234678-HpCDD	29.747	0.4449	0.112	0.01	0.001	Valid
1234789-HpCDF	30.475	0.0453	0.096	0.01	0.000	nd<LOQc
OCDD	34.197	0.4003	0.112	0.0003	0.000	Valid
OCDF	34.457	0.0417	0.096	0.0003	0.000	nd<LOQc
Sum					1.163	pg/g


Fish sample WHO-TEQ (2005) – 78.6g Fresh weight – HR-GCMS 1.29pg/g

Sample Chromatography And Calibration – PCDD/Fs



Fish sample WHO-TEQ (2005) 1.16pg/g – 78.6g Fresh weight

Standard Chromatography And Calibration – dl-PCBs (WP-CS1)



Dioxin Analysis Legislation - Approaches to Dioxin Analysis

- Europe
 - Performance based approach.
 - Regulated technology and defined performance and quality criteria.
 - Methods are provided as a guide and can be modified.
 - Levels in food are regulated with strict maximum limits.
 - Emission levels are regulated.
- North America & Japan
 - Method based approach
 - Regulatory agencies (such as USEPA) provide complete methods to be followed.
 - Deviations to the method are possible. But must be documented closely.
 - Significant deviations from methods can result in the analysis being challenged or becoming unsellable in market.
 - Levels in food are not regulated.
 - Emission guidelines are given.
- ROW adopt variety of approaches,
- dependent on regulatory and/or business requirement

- **Confirmatory Analysis**

- Unequivocal identification and quantification
- Results can be used in prosecution



- **Screening Analysis**

- Indication of levels only
- Must be followed with confirmatory analysis



What Has Been Happening In EU Dioxin Analysis?

ANALYTICAL CRITERIA FOR USE OF MS/MS FOR DETERMINATION OF DIOXINS AND DIOXIN-LIKE PCBS IN FEED AND FOOD

Kotz A¹, Malisch R^{1*}, Focant J², Eppe G², Cederberg TL³, Rantakokko P⁴, Fürst P⁵, Bernsmann T⁵, Leondiadis L⁶, Lovász C⁷, Scortichini G⁸, Diletti G⁸, di Domenico A⁹, Ingelido AM⁹, Traag W¹⁰, Smith F¹¹, Fernandes A¹¹

¹European Union Reference Laboratory (EU-RL) for Dioxins and PCBs in Feed and Food, Freiburg, Germany; ²CART, University of Liege, Belgium; ³National Food Institute, Technical University of Denmark, Søborg, Denmark; ⁴National Institute for Health and Welfare (THL), Kuopio, Finland; ⁵Chemisches und Veterinäruntersuchungsamt Münsterland-Emscher-Lippe (CVUA-MEL), Münster, Germany; ⁶NCSR Demokritos, Athens, Greece; ⁷Central Agricultural Office, Budapest, Hungary; ⁸Istituto Zooprofilattico Sperimentale dell'Abruzzo e del Molise "G. Caporale, Teramo, Italy; ⁹Istituto Superiore di Sanità (ISS), Roma, Italy; ¹⁰RIKILT, Wageningen, The Netherlands; ¹¹The Food and Environment Research Agency (FERA), York, United Kingdom

<http://www.dioxin20xx.org/pdfs/2012/1041.pdf>

- 2009–2012 GC-MS/MS systems were evaluated from several vendors
- Conclusion:
 - “GC-MS/MS systems with sufficient sensitivity” can be “used as confirmatory methods”
 - “However, for the determination of low background levels the application of GC-HRMS is required” e.g. <20% of Maximum Level (ML)

Does MS/MS Agree with HRMS?




Figure 11. Comparative results (upperbound concentration values) for 50 food and feed samples analyzed by GC/HRMS and GC/MS/MS.


- At EU ML results are comparable.
- Below EU ML results are not comparable.

Does MS/MS Agree with HRMS?

Comparison of GC-MS/MS Results of laboratories (3)

PCDD/Fs

- Chemisches und Veterinäruntersuchungsamt (CVUA) Freiburg (EU-RL for Dioxins and PCBs), GC-MS/MS



Comparable results for MS/MS and HRMS for WHO-TEQ, deviation of MS/MS mostly below 20 %

14 / GC-MS/MS

EURL
European Union Reference Laboratory
Dioxins and PCBs

State Institute for Chemical and
Veterinary Analysis of Food
CVUA Freiburg



- At EU ML results are comparable.
- Below EU ML results are not comparable.

What's New In EU Regulation For Dioxins?

Past

[Commission Regulation \(EU\) No 252/2012 of 21 March 2012](#)

- Specifies use of **GC-HRMS** for confirmatory dioxin analysis
- **GC-MS/MS** was allowed a screening technique.



Present

[Commission Regulation \(EU\) No 589/2014 of 2 June 2014](#)

- Specifies use of **GC-HRMS** or **GC-MS/MS** for confirmatory dioxin analysis
- **GC-MS/MS** is “an appropriate confirmatory method for checking compliance with the maximum”, only.
- **GC-HRMS** remains the recommended technique for “determination of low background levels in food monitoring, following of time trends, exposure assessment of the population”.

Compliance With New EU Dioxin Regulation

European Commission GC-MS/MS Confirmatory Performance Criteria	TSQ 8000 Evo GC-MS/MS Capabilities	Compliance Confirmed
Two specific precursor ions with two specific production ions	All recommended methods developed as defined in criteria Fully optimized by AutoSRM	Yes
Tolerance of ion ratios within $\pm 15\%$	< 10% measured at EPA 1613 CSL level (n=14)	Yes
Resolution of each quadrupole equal to or better than unit mass resolution	All recommended methods developed Q1 and Q3 at 0.7 Da	Yes
The % RSD of the five (or more) Relative Response Factors (RRFs) for each unlabeled PCDD/PCDF and labelled internal standards must not exceed 20%	6 point curve EPA 1613 CSL-CS4 2%	Yes

Applications and Methods

Application		Norm	Thermo Scientific TSQ 8000 Evo GC-MS/MS	Thermo Scientific DFS GC-HRMS
Food	EU Regulation	EU Regulatory Feed Control (at ML)	Approved	Approved
Food	EU Methods	EN 16215	Not in method	Approved
Food	EURL Recommendation	Background food studies (<1/5th EU ML)	Not recommended	Recommended
Clinical	EURL Recommendation	Human studies at trace levels	Not recommended	Recommended
Environmental	US Method	US EPA 8280A	Approved	Not in method
Environmental	EU Method	EN 1948	Not in method	Approved
Environmental	US Method	US EPA 1613 B for strict EPA compliance	Not in method	Approved
Environmental	US Method	US EPA Method 23	Not in method	Approved
Environmental	US Method	US EPA Method 8290	Not in method	Approved
Environmental	Japan Method	JIS K0311	Not in method	Approved
Environmental	Japan Method	JIS K0312	Not in method	Approved

- Any laboratory following EU approach will now be able to use GC-MS/MS for ML compliance control
- Laboratories strictly following US/Japan methods and regulation will continue to use GC-HRMS