


Primer selection guidelines utilizing VeriQuest® qPCR for microarray validation

This tutorial provides examples of appropriate primer design strategies when using USB® VeriQuest qPCR master mixes to validate Affymetrix® microarrays. These guidelines generate optimal concordance of microarray and qPCR results when using **First-Strand cDNA Synthesis Kit (PN 75780)** followed by real-time PCR with **VeriQuest Probe (PN 75650)** or **VeriQuest SYBR® Green (PN 75600)** qPCR master mixes.

Note: Assays can be performed using either VeriQuest SYBR Green or Probe chemistries. SYBR assays are presented here to describe the most cost-effective solution.

Example 1: Validating GeneChip® Human Transcriptome Array (HTA) 2.0 data


GeneChip HTA 2.0 is the highest resolution microarray for whole transcriptome expression profiling of all transcript isoforms. A key advantage of HTA 2.0 is a simple, fast, and free data analysis package consisting of **Expression Console™ Software** and **Transcriptome Analysis Console (TAC) Software** that can be downloaded from the Affymetrix website (Affymetrix.com). Expression Console generated .chp files can be opened in TAC for analysis and easy interpretation of differential gene-level expression as well as alternative splicing between samples assayed on HTA 2.0.




The following workflow demonstrates how to design VeriQuest qPCR assays for validating alternative splicing. HTA 2.0 data for Human Total RNA from liver and brain were subjected to Alternative Splicing Analysis in TAC to obtain a list of differentially expressed genes and alt-splicing events (liver vs. brain). To increase the confidence of the TAC analysis, a p value of 0.001 was used, which is more stringent than the default value of 0.05. Significant alt-splicing events were selected based on **Probe Selection Region (PSR)**

Splicing Index (SI) values. In this example, the Clathrin Light Chain B (CLTB) gene, which helps regulate receptor-mediated endocytosis, was selected based on its criteria of having a low **Gene Fold Change** = <2 or >2 (**-1.37**) and a significant **PSR SI** value = >2 or <-2 (**-6.11**). The **Junction Probe (JUC) SI** values (**-8.89**, **-5.64**) provide further supporting evidence for differential splicing.

By clicking on the row containing CLTB in TAC, all known transcript isoforms in the database are displayed. In this gene structure view, significant alt-splicing events are darker in color and will highlight when clicking on the respective PSR as shown below (view filtered to show only the accession numbers for isoforms in RefSeq). In this example, CLTB exon 8 is characterized as a Cassette Exon that is skipped more often in liver versus brain.



To validate the CLTB cassette exon alt-splicing event with VeriQuest qPCR, design multiple primer sets (2-4) within constitutive exons (exons that are expressed in both conditions or always “on”) and multiple primer sets (2-4) within the alt-spliced exon. In this example, CLTB exons 1, 2, 6, 7, and 9 are constitutive and exon 8 is alternatively spliced. If the exon length is less than 70 bp, as for CLTB exon 8, design qPCR primers with one of the pair within the alt-spliced exon and the other within the flanking constitutive exon as shown below. Real-time PCR data from constitutive (C) amplicons provide gene-level expression information. Data from alt-spliced (A) amplicons are used to calculate qPCR Splicing Indexes. This strategy enables accurate calculation of Splicing Index from VeriQuest qPCR data (**Appendix A**) for optimal array concordance. The figure below shows the CLTB gene structure view from TAC (to scale) and approximate positions of qPCR amplicons (not to scale) designed to validate HTA 2.0 data.



Designing primers for exons greater than 70 bp (optimal amplicon length 70-100 bp)

Designing qPCR primers to interrogate specific exons is simplified with TAC. By right clicking on an exon of interest (CLTB constitutive exon 7 below), and selecting **Copy DNA Sequence (5' to 3')**, the sequence can be copied and pasted into **NCBI Primer-BLAST** for efficient design of qPCR primers for use with VeriQuest SYBR Green qPCR Master Mix.

The screenshot shows the TAC interface with a genomic view of the CLTB gene. A red circle highlights a context menu that appears when right-clicking on exon 7. The menu options are: Copy Position, Copy DNA Sequence (5' to 3'), View Exon in UCSC, and View Exon in ENSEMBL.

NCBI can be opened from TAC by right clicking on the RefSeq accession number and selecting **Link to NCBI**:


The screenshot shows the TAC interface with a genomic view of the CLTB gene. A red circle highlights a context menu that appears when right-clicking on the RefSeq accession number NM_001834. The menu options are: Copy ID, Link to Ensembl, and Link to NCBI.

Paste the exon sequence into **NCBI Primer-BLAST**, set Max PCR product size to 100, and click **Get Primers**:

The screenshot shows the NCBI Primer-BLAST interface. The PCR Template section shows the sequence for CLTB 5'-3'. The Primer Parameters section shows the Max PCR product size set to 100 and the number of primers to return set to 10.

Designing primers for exons less than 70 bp (optimal amplicon length 70-100 bp)

To interrogate exons less than 70 bp, as is the case for CLTB alt-spliced exon 8 (54 bp), design qPCR primers spanning exon-exon junctions by right clicking on the RefSeq accession number that contains the alt-spliced exon 8 (in this case NM_007097) and selecting **Link to NCBI**:



Click the **Nucleotide** link under Genomes:

Genomes			Chemicals		
Assembly	0	genomic assembly information	BioSystems	0	molecular pathways with links to genes, proteins and chemicals
BioProject	0	biological projects providing data to NCBI	PubChem BioAssay	0	bioactivity screening studies
BioSample	0	descriptions of biological source materials	PubChem Compound	0	chemical information with structures, information and links
Clone	0	genomic and cDNA clones	PubChem Substance	0	deposited substance and chemical information
dbVar	0	genome structural variation studies			
Epigenomics	0	epigenomic studies and display tools			
Genome	0	genome sequencing projects by organism			
GSS	0	genome survey sequences			
Nucleotide	1	DNA and RNA sequences			
Probe	20	sequence-based probes and primers			
SNP	0	short genetic variations			

This will open the NCBI Reference Sequence for CLTB:

NCBI Resources How To Sign in to NCBI

Nucleotide Nucleotide Search Advanced Help

Display Settings: GenBank Send: Change region shown Customize view

Homo sapiens clathrin, light chain B (CLTB), transcript variant 2, mRNA

NCBI Reference Sequence: NM_007097.3

FASTA Graphics

Go to: Analyze this sequence Run BLAST Pick Primers

LOCUS NM_007097 1280 bp mRNA linear PRI 15-MAR-2015

The exon sequences can be viewed towards the bottom of the page and start/end positions determined.

exon

exon

STS

STS

STS

STS

polyA site

/inference="alignment:Splice:1.39.8"

670..723

/gene="CLTB"

/gene_synonym="LCB"

/inference="alignment:Splice:1.39.8"

724..1214

/gene="CLTB"

/gene_synonym="LCB"

/inference="alignment:Splice:1.39.8"

872..1111

/gene="CLTB"

/gene_synonym="LCB"

/standard_name="RH98811"

/db_xref="UniSTS:84078"

896..1170

/gene="CLTB"

/gene_synonym="LCB"

/standard_name="WI-7674"

/db_xref="UniSTS:58080"

968..1161

/gene="CLTB"

/gene_synonym="LCB"

/standard_name="RH69093"

/db_xref="UniSTS:58244"

987..1133

/gene="CLTB"

/gene_synonym="LCB"

/standard_name="G54012"

/db_xref="UniSTS:109394"

1214

/gene="CLTB"

/gene_synonym="LCB"

ORIGIN

1

ccgctctctcc

cgcgcgcggcg

ccgcacccgc

agtgacagcc

agcggggccc

ggtggcggag

61

aggaagtgcg

gtccgcgccca

agcccgctcc

cgcgcagccc

ggctccccgc

ggctcgggtg

121

acagcgtcgc

ggccgcgcgga

cgcagcgcgg

ggcaggcgcg

ggcagagccg

agcgcagcgg

181

aggctccggc

ggaggcgccg

ggaaaaatgg

tgatgacttt

ggctttctct

cgctcgtcga

241

gagcgggtgc

cgggagcgcg

cggagggagg

cccgggcgcc

gccttctctg

cccagcagga

301

gagcgagatt

gcaggcatag

agaacgcaga

gggcttcggg

gcacctgcgc

gcagccatgc

361

ggcccccgcg

cagccggggcc

ccacgagtgg

ggctgggtct

gaggacatgg

ggaccacagt

421

caatggagat

gtgtttcagg

aggccaacgg

tcctgctgat

ggctacgcag

ccattgcccc

481

ggctgacagg

ctgaccacag

agcctgagag

catccgcaag

tggcgagagg

agcagaggaa

541

acggctgcaa

gagctggatg

ctgcattctaa

ggtcacggaa

caggaatggc

gggagaaggg

601

caagaaggac

ctggaggagt

ggaaccacgc

ccagagttaa

caagtagaga

agaacaagat

661

caacaaccgg

atcgctgaca

aagcattcta

ccagcagcca

gatgctgata

tcacgggcta

721

cgtggcatcc

gaggaggctt

tcgtgaagga

atccaaggag

gagaccccag

gcacagagtg

781

ggagaagggt

gcccagctat

gtgacttcaa

ccccaaagag

agcaagcagt

gcacaagatg

841

gtcccccgct

cgctcgggtc

tcattgtccc

gaagcagacg

ccactgtccc

gctaggtgcc

901

tgctagggtc

atggccacag

agcatgggct

gggcctgggc

acaggaggag

cagctgcttt

961

ggtcgggggt

gagactcgca

gcagctgcta

cccacagcct

attccactcc

tccccatctc

1021

caggcgctgg

gagggggggc

ctcaccccat

cacgcctcgc

tcctctctgg

ccctctgggt

1081

cagccctcca

cgctcctctc

cagtctactc

aattgtgact

gtccctctcg

atgtattttt

1141

tttcttgaat

taaaagatgt

atttttgaat

ctttttacac

ttattttatta

tcatttctac

CLTB NM_007097

exon	start	end	length
1	1	392	392
2	393	439	47
6	440	557	118
7	558	669	112
8*	670	723	54
9	724	1214	491

* CLTB exon 8 alt-spliced as determined by TAC

CLTB exon 8 sequence selected in NCBI. This sequence can be cross-referenced with TAC for confirmation of the correct alt-spliced exon.

By clicking **Pick Primers** under Analyze this sequence, qPCR primer sets can be designed using **NCBI Primer-BLAST**.

NCBI Resources How To Sign in to NCBI

Nucleotide Nucleotide Search Advanced Help

Display Settings: GenBank Send: Change region shown Customize view

Homo sapiens clathrin, light chain B (CLTB), transcript variant 2, mRNA

NCBI Reference Sequence: NM_007097.3

FASTA Graphics

Go to:

LOCUS

NM_007097

1280 bp

mRNA

linear

FRI 15-MAR-2015

DEFINITION

Homo sapiens clathrin, light chain B (CLTB), transcript variant 2, mRNA.

ACCESSION

NM_007097

VERSION

NM_007097.3

GI:365906278

KEYWORDS

RefSeq.

SOURCE

Homo sapiens (human)

ORGANISM

Homo sapiens

Analyze this sequence

Run BLAST

Pick Primers

Highlight Sequence Features

Find in this Sequence

Articles about the CLTB gene

Genome-wide cDNA screening identifies hnet

5

To design qPCR primers that span exon-exon junctions with one of the pair within the exon of interest (in this case CLTB alt-spliced exon 8), enter the start/end positions of exon 7 (558/669) in the **Forward primer From/To** fields and the start/end positions of exon 8 (670/723) in the **Reverse primer From/To** fields. Set Max PCR product size to 100, and click **Get Primers**.

Primer-BLAST: A tool for finding specific primers

NCBI/ Primer-BLAST: Finding primers specific to your PCR template (using Primer3 and BLAST).

Reset page Save search parameters Retrieve recent results Publication Tips for finding specific primers

PCR Template

Enter accession, gi, or FASTA sequence (A refseq record is preferred) [Clear](#)

NM_007097.3

Or, upload FASTA file [Browse...](#)

Range

Forward primer From To

Reverse primer From To [Clear](#)

558 669

670 723

Primer Parameters

Use my own forward primer (5'→3' on plus strand) [Clear](#)

Use my own reverse primer (5'→3' on minus strand) [Clear](#)

PCR product size

Min Max

70 100

of primers to return


10

Primer melting temperatures (T_m)

Min Opt Max Max T_m difference

57.0 60.0 63.0 3 [Clear](#)

Select a qPCR primer pair (optimal amplicon length 70-100 bp) for interrogating the CLTB alt-spliced exon 8.



Checking primers for specificity and SNPs

Primer pairs can be checked for specificity in **NCBI Primer-BLAST** by inputting the forward and reverse sequences in the **Use my own primer** fields and clicking **Get Primers**:

Primer-BLAST: A tool for finding specific primers

NCBI/ Primer-BLAST: Finding primers specific to your PCR template (using Primer3 and BLAST).

Reset page Save search parameters Retrieve recent results Publication Tips for finding specific primers

PCR Template

Enter accession, gi, or FASTA sequence (A refseq record is preferred) [Clear](#)

NM_007097.3

Or, upload FASTA file [Choose File](#) No file chosen

Range

Forward primer From To

Reverse primer From To [Clear](#)

Primer Parameters

Use my own forward primer (5'→3' on plus strand) [Clear](#)

Use my own reverse primer (5'→3' on minus strand) [Clear](#)


PCR product size

Min Max


AGCGCCAGAGTGAACAAGTAG

CCGATGATATCAGCATCTGGCT

qPCR primer pairs that are unique for a specific target should only give products on the target of interest (in this case CLTB exon 7/8) of the appropriate amplicon length (in this case 91).



UCSC In-Silico PCR can also be used to check qPCR primer pairs for specificity and SNP positions. Common SNPs within primer annealing sites should be avoided. The UCSC database can be opened from TAC by right clicking on any exon and selecting **View Exon** in UCSC.



When primer pairs span exon-exon junctions, a longer amplicon may be amplified due to genomic DNA contamination and will show as a higher molecular weight band on agarose gels. This contamination will interfere with SYBR-based qPCR analysis and should be minimized by using isolation kits with a DNase step (e.g. [PrepEase RNA Spin Kit, PN 78766](#)) or by DNase treatment of isolated RNA using [rDNase I, RNase-Free \(PN 78411\)](#). Genomic DNA contamination is more difficult to detect when primer pairs are within an exon boundary, and the inclusion of a no reverse transcriptase control is recommended in those situations.

Fold change from real-time PCR data is calculated using the $\Delta\Delta C_t$ method relative to an endogenous reference gene (e.g. ACTB). Gene-level fold change is calculated using the fold change of constitutive exons. Splicing Index from qPCR data is calculated for alt-spliced amplicons by normalizing their fold change to the average fold change of constitutive amplicons. In effect, Splicing Index is a fold-change measure of the alt-splice event. See **Appendix A** for detailed calculation and table of recommended reference genes and primers in **Appendix B**.


Validation of Affymetrix GeneChip HTA 2.0 alternative splicing data by VeriQuest SYBR Green qPCR assays may have differences in the magnitude of the splicing indexes. This is normal due to the inherent differences in the dynamic range and noise of each platform. In an internal study conducted by Affymetrix scientists, 16 separate differential alt-splicing events selected from HTA 2.0 data were confirmed by VeriQuest qPCR with 100%-directional concordance (see [White Paper](#)).

Example 2: Validating GeneChip Human Gene (HuGene) 2.0 ST Array

When using VeriQuest qPCR Master Mixes to validate Affymetrix 3' IVT (e.g. Human Genome U133 Plus 2.0) and Gene (e.g. Human, Mouse, Rat Gene 1.x or 2.x ST) microarrays, qPCR primers should be designed against constitutive exons using a workflow similar to **Example 1**. While these arrays are not annotated to the extent of HTA 2.0 in **TAC, Expression Console** generated files can be analyzed in TAC Software for easy interpretation of differential expression between samples.

In this example, HuGene 2.0 data from Human liver and muscle total RNA were analyzed by TAC Software to obtain a list of differentially expressed genes (see figure below). The Complement Factor H-Related 3 (CFHR3) gene, which encodes a protein that binds heparin and may be involved in complement regulation, was selected based on its significant gene-level fold change between samples (liver vs. muscle) and constitutive regions of significant length for interrogation by VeriQuest qPCR. The highlighted CFHR3 gene's differential expression is shown as a blue 'x' in the displayed scatter plot (higher expression in liver vs. muscle). The **Transcript Cluster ID** (16675428) can be queried in **NetAffx™ Analysis Center**, which is accessible on the Affymetrix website ([Affymetrix.com](#)), for appropriate primer design strategy.


Right click the Transcript Cluster ID for CFHR3 and select **Search Affymetrix NetAffx**:



Click the **Entrez ID (10878)** under Transcripts Detected By 16675428 to link out to **NCBI** for the gene structure view:

Transcripts Detected By 16675428				
Full Length	Gene Symbol (Entrez ID)	Description	Assign Score	Coverage
		Pathways		
BC070259 GenBank	CFHR3 (10878)	Homo sapiens complement factor H-related 3, mRNA (cDNA clone IMAGE:4716670), complete cds.	100 (8/8)	67 (8/12)
ESTs	Gene Symbol (Entrez ID)	Description	Assign Score	Coverage
		Pathways		
Predicted	Gene Symbol (Entrez ID)	Description	Assign Score	Coverage
		Pathways		


From the NCBI gene structure view of CFHR3, constitutive exons (e.g. CFHR3 exons 1-3, 5,6) are present in all variants where potentially alt-spliced exons (e.g. CFHR3 exon 4) are not. Use constitutive exon sequences for designing qPCR primers for gene-level validation in VeriQuest qPCR assays.



Note: Gene structures are subject to change with discovery of new variants. Non-concordance between assays due to undefined alt-splicing events can be resolved with a 2nd qPCR assay with primer pairs in a different constitutive region.

The RefSeq accession number (NM_021023) can then be queried on the [NCBI](#) website to determine sequences of constitutive exons for qPCR primer design (similar to the workflow in **Example 1**) for validation with VeriQuest SYBR Green qPCR Master Mixes.

As an alternative to primer design using [NCBI Primer-BLAST](#) and [UCSC In-Silico PCR](#), both SYBR- and Probe-based pre-designed qPCR assays from [IDT](#) can be used for gene-level validation of Affymetrix expression arrays. By simply inputting the **Gene Symbol**, **Species**, and **Assay Configuration** preference (in this case CFHR3, Human, Intercalating dyes, primers only) and clicking **Submit**, available assays for **Detecting All Variants** can be purchased. Notice there are no assays available for exon 4 which is a non-constitutive exon and should be avoided in VeriQuest qPCR validation assays.



Sign In

Chat Now

0 Items \$0.00

Order Menu

Products & Services

Support & Education

Tools

Search

Go

PrimeTime[®] qPCR Assays

IDT now offers PrimeTime qPCR Assays that are guaranteed to work for human, mouse, and rat transcriptomes. All PrimeTime qPCR Assays consist of two primers and a hydrolysis probe. All three components are combined into a single tube and shipped in 2–4 business days. Each oligo undergoes 100% QC by mass spectrometry, with all QC results provided free of charge on the IDT website.

Primers Alone. Primers for these assays are also available without probe, for use with intercalation dyes such as SYBR[®] Green and EvaGreen[®]. For primers alone, under **Default Assay Configuration**, choose **Intercalating dyes, Primers only**.

All assays are designed to detect gene expression. To design assays for copy number variation, please [contact us](#).

Basic

Batch

Help

Gene Symbol

RefSeq

Assay ID

Cfhr3

+

+

+

Species

Human

Mouse

Rat

Exact Match

Default Assay Configuration

5' nuclease, probe included

Intercalating dyes, primers only

Search

Add to Order

Displaying 5 results for > Gene: Cfhr3

Export

Filter

Compare

Assay ID	Gene Query	NCBI Gene Symbol	Ref Seq #	Detects All Variants	Exon Location	Assay Configuration
Hs.PT.58.1202575	CFHR3	CFHR3	NM_001166624(2)	Yes	3 - 5	Std, DNA Primer
Hs.PT.58.25127850.g	CFHR3	CFHR3	NM_021023(2)	Yes	2 - 2	Std, DNA Primer
Hs.PT.58.25506640.g	CFHR3	CFHR3	NM_021023(2)	Yes	3 - 3	Std, DNA Primer
Hs.PT.58.3291139.g	CFHR3	CFHR3	NM_001166624(2)	Yes	2 - 3	Std, DNA Primer
Hs.PT.58.40280646.g	CFHR3	CFHR3	NM_001166624(2)	Yes	6 - 6	Std, DNA Primer

Fold change from real-time PCR data is calculated using the $\Delta\Delta C_t$ method relative to an endogenous reference gene (e.g. ACTB). Gene-level fold change is the fold change of constitutive exons.

Appendix A

The Splicing Index for HTA data represents the ratio of the exon intensities in condition 1 versus condition 2 after normalization to their respective gene intensities in each sample. A Splicing Index value of 0 indicates that the Probe Selection Region (PSR) is present at equal levels in both conditions. A positive Splicing Index value implies elevated inclusion, and a negative value suggests increased PSR skipping in condition 1 versus condition 2.

This is an example of fold change and Splicing Index calculation from VeriQuest qPCR data. See **Example 1** above for CLTB gene structure view and primer design.

														Fold Change		Splicing Index	
Gene	Amplicon		C _t per triplicates						Average C _t		Average ΔC _t (ACTB-CLTB)		ΔΔC _t	2 ^{ΔΔC_t}	2 ^{ΔΔC_t} or -1/2 ^{ΔΔC_t}	FC(A)/AvgFC(C)	SI
			Liver			Brain			L	B	L	B					
CLTB	constitutive (C)	exon 1-1	31.3	31.2	31.3	29.9	29.8	29.8	31.3	29.8	-4.96	-3.95	-1.01	0.50	-2.0		
		exon 7-7	29.9	29.8	30.0	28.8	28.9	28.9	29.9	28.9	-3.56	-3.00	-0.56	0.68	-1.5		
		exon 9-9	31.5	31.2	31.3	29.9	29.9	29.9	31.3	29.9	-5.03	-4.02	-1.01	0.50	-2.0		
	alt-spliced (A)	exon 7-8	35.0	35.0	35.0	29.1	29.0	29.1	35.0	29.0	-8.70	-3.14	-5.55	0.021		0.038	-26
		exon 8-9	35.0	35.0	35.0	28.9	28.9	28.9	35.0	28.9	-8.70	-2.99	-5.70	0.019		0.034	-29
ACTB	constitutive	exon 3-4	26.3	26.3	26.3	25.9	25.9	25.9	26.3	25.9							

Fold change and Splicing Index from qPCR data are calculated using the $\Delta\Delta C_t$ method relative to an endogenous reference gene (ACTB in this example). A raw C_t of 35 is used as the limit of detection: C_t values are set at 35 for any replicates with C_t values not determined or >35.

Calculate average ΔC_t as the average C_t (ACTB) - average C_t (CLTB) for both conditions (liver and brain).

Calculate $\Delta\Delta C_t$ as ΔC_t (liver) - ΔC_t (brain).

Calculate 2 $^{\Delta\Delta C_t}$.

Gene-level fold change is reported as the fold change of constitutive exons (C).

Fold change = 2 $^{\Delta\Delta C_t}$ if 2 $^{\Delta\Delta C_t} \geq 1$

Fold change = -1/2 $^{\Delta\Delta C_t}$ if 2 $^{\Delta\Delta C_t} < 1$

For the constitutive exon 1, the calculated 2 $^{\Delta\Delta C_t}$ value is less than 1 (0.50), indicating slightly decreased gene expression in liver versus brain. This is finally reported as -1/0.50 = -2.0, as a negative number.

Splicing Index from qPCR data is calculated for alt-spliced amplicons (A) by normalizing their fold change to the average fold change of constitutive exons.

Calculate FC (A)/average FC (C)

For amplicon spanning exons 7-8:

0.021/average (0.50, 0.68, 0.50) = 0.038

Splicing Index = FC (A)/average FC (C) if FC (A)/average FC (C) ≥ 1

Splicing Index = -1/(FC (A)/average FC (C)) if FC (A)/average FC (C) < 1

For amplicon spanning exons 7-8, the calculated FC (A)/average FC (C) value is less than 1 (0.038), indicating decreased exon 8 inclusion in liver. This is finally reported as -1/0.038 = -26, as a negative number.

Appendix B

Real-time PCR data is normalized to a reference gene, an mRNA that is stably expressed under the experimental conditions employed. It is recommended to screen multiple reference genes to determine which expression levels fluctuate the least. Commonly used endogenous reference genes with possible forward and reverse primer sequences are listed in the table.

Gene ID	Description	Primer sequences (5' to 3')	
18S	18S ribosomal RNA	F	CGAAGACGATCAGATACCGT
		R	GGTCATGGGAATAACGCCG
28S	28S ribosomal RNA	F	TCGTCCGACCTGGGTATAG
		R	GCTATCCTGAGGGAACTTCG
ACTB	actin, beta	F	AACCGCGAGAAGATGACCCAGAT
		R	TAGCACAGCCTGGATAGCAACGTA
GAPDH	glyceraldehyde-3-phosphate dehydrogenase	F	ATCGTGGAAGGACTCATGACCACA
		R	TAGAGGCAGGGATGATGTTCTGGA
HPRT1	hypoxanthine phosphoribosyltransferase 1	F	ATGGACAGGACTGAACGTCTTGCT
		R	TTGAGCACACAGAGGGCTACAATG
POLR2A	polymerase (RNA) II (DNA directed) polypeptide A	F	GCATTGACTTGCGTTTCCAC
		R	CCATCACACATGTGCCGTTC
PPIA	peptidylprolyl isomerase A (cyclophilin A)	F	GTCTCCTTTGAGCTGTTTGC
		R	AAGCAGGAACCCCTTATAACC
TBP	TATA box binding protein	F	ACCAGGTGATGCCCTTCTGTAAGT
		R	ATGAGCAACTCACAGTCACGCT

For Research Use Only. Not for use in diagnostic procedures.

Affymetrix, USB, GeneChip, and VeriQuest are registered trademarks of Affymetrix, Inc. NetAffx is a trademark of Affymetrix, Inc. SYBR is a registered trademark of Molecular Probes, Inc. and is provided under an agreement with Molecular Probes, Inc. All trademarks are property of their respective owners.

Affymetrix, Inc.

USB® Products

26111 Miles Road
Cleveland, Ohio 44128
Tel: 800-321-9322 | 216-765-5000
Fax: 800-535-0898 | 216-464-5075
USBcustomerserv@affymetrix.com
usb.affymetrix.com

Affymetrix UK Ltd.

USB® Products

Voyager, Mercury Park,
Wycombe Lane, Wooburn Green,
High Wycombe HP10 0HH, UK
Tel: +44 (0)1628 55 2600
Fax: +44 (0)1628 55 2675
USBcustomerserveurope@affymetrix.com

Affymetrix Pte Ltd.

USB® Products

7 Gul Circle, #2M-01
Keppel Logistics Building
Singapore 629563
Tel: +65 63957310
Fax: +65 63957300
USBcustomerserv@affymetrix.com