


# High Sensitivity ELISA

## For Human and Mouse



# High Sensitivity ELISA Kits



## Low Abundance Target Analysis Made Easy


High Sensitivity ELISAs from eBioscience enable the analysis of target analytes in low concentrations from serum, plasma, or cell culture samples. Scientists who are looking for improved sensitivity can now take advantage of our High Sensitivity ELISA Portfolio.

eBioscience High Sensitivity ELISA Kits offer accurate quantification below conventional ELISA limits. This is made possible by the inherent design which amplifies the signal generated by the enzyme horseradish peroxidase (HRP) and the colorimetric-detectable substrate tetramethylbenzidine (TMB), based on Tyramide Signal Amplification (TSA) technology\*.

\*High Sensitivity ELISAs are based on a patent from PerkinElmer LAS, patent number 5,731,158; 5,583,001; 5,196,306.

## Advantages of High Sensitivity ELISAs

- **Sensitive** – Detect low cytokine concentrations  $\leq 1.0$  pg/mL for reliable quantification.
- **Low Sample Volume** – 50  $\mu$ L sample volume required. Ideal for limited sample resources.
- **Ease-of-Use** – Straightforward protocol with all required reagents included.
- **Existing Instrumentation** – Uses standard colorimetric plate reader (absorption measurement at 450 nm).
- **Includes reagents** for 96-well microtiter plate.







Human				
Analyte	Sensitivity	Standard Range	Cat. No.	Sample Volume
IFN- $\gamma$	0.06 pg/mL	0.16 - 10.0 pg/mL	BMS228HS	50 $\mu$ L
IL-1 $\beta$	0.05 pg/mL	0.16 - 10.0 pg/mL	BMS224HS	50 $\mu$ L
IL-2	0.4 pg/mL	0.9 - 60.0 pg/mL	BMS221HS	50 $\mu$ L
IL-4	0.1 pg/mL	0.25 - 16.0 pg/mL	BMS225HS	50 $\mu$ L
IL-6	0.03 pg/mL	0.08 - 5.0 pg/mL	BMS213HS	50 $\mu$ L
IL-10	0.05 pg/mL	0.39 - 25.0 pg/mL	BMS215HS	50 $\mu$ L
IL-12 (p70)	0.1 pg/mL	0.16 - 10 pg/mL	BMS238HS	50 $\mu$ L
IL-17A	0.1 pg/mL	0.23 - 15 pg/mL	BMS2017HS	50 $\mu$ L
TNF- $\alpha$	0.13 pg/mL	0.31 - 20.0 pg/mL	BMS223HS	50 $\mu$ L

Mouse				
Analyte	Sensitivity	Standard Range	Cat. No.	Sample Volume
IFN- $\gamma$	0.05 pg/mL	0.31 - 20 pg/mL	BMS606HS	50 $\mu$ L
IL-2	0.05 pg/mL	0.78 - 50 pg/mL	BMS601HS	50 $\mu$ L
IL-4	0.32 pg/mL	0.31 - 20 pg/mL	BMS613HS	50 $\mu$ L
IL-6	0.21 pg/mL	1.56 - 100 pg/mL	BMS603HS	50 $\mu$ L

## Look Beyond the Borders of Conventional Detection Limits

The eBioscience High Sensitivity ELISA provides a solution for difficult to detect, low expressing proteins by further enhancing the signal obtained from conventional colorimetric ELISA. Whereas conventional sandwich ELISA's amplification is accomplished only via the biotin/streptavidin-HRP amplification system, the eBioscience High Sensitivity ELISAs further enhance this signal with the additional signal amplification step using a Tyramide amplification reagent. The biotinylated Tyramide (Amplification Reagent I) undergoes a conversion into a short-lived, extremely reactive intermediate in the presence of HRP. These free radicals bind covalently to any adjacent protein in the vicinity of HRP. This generates multiple biotin binding sites for additional streptavidin-HRP molecules.

- The amount of reacted biotinyl-Tyramide is proportional to the amount of HRP in the well.
- Following incubation, unbound biotinylated Tyramide is removed during a wash step.
- After removal of excess biotin-Tyramide from the wells, streptavidin-HRP (Amplification Reagent II) is added, and binds to the additional biotin sites created during the biotin-Tyramide reaction. This amplification step effectively multiplies the quantity of HRP molecules for the final TMB substrate reaction.



	Standard ELISA		High Sensitivity ELISA	
	Range pg/mL	Sensitivity pg/mL	Range pg/mL	Sensitivity pg/mL
Human TNF-α	8 - 500	2.3	0.31 - 20	0.13
Human IFN-γ	1.6 - 100	0.99	0.16 - 10	0.06
Human IL-1β	3.9 - 250	0.32	0.16 - 10	0.05
Human IL-2	19 - 1,200	9.9	0.94 - 60	0.40
Human IL-4	8 - 500	1.3	0.25 - 16	0.10
Human IL-6	1.6 - 100	0.9	0.08 - 5	0.03
Human IL-10	3.2 - 200	0.99	0.39 - 25	0.05
Human IL-12p70	3.2 - 200	2.1	0.16 - 10	0.10
Human IL-17A	1.6 - 100	0.5	0.23 - 15	0.10
Mouse IL-2	15.6 - 1,000	5.3	0.78 - 50	0.05
Mouse IL-4	3.9 - 250	1.98	0.31 - 20	0.32
Mouse IL-6	31.3 - 2,000	6.5	1.56 - 100	0.21
Mouse IFN-γ	15.6 - 2,000	5.3	0.31 - 20	0.05

**Headquarters**

eBioscience, Inc.  
10255 Science Center Drive  
San Diego, CA 92121  
USA

**SERVICE AND SUPPORT FOR DIRECT SALES**

**North America**

Technical Support:  
*For Research Products:*  
888.810.6168  
858.642.2058  
tech@eBioscience.com  
*For Clinical Products:*  
877.726.8559  
858.642.2058  
tech@eBioscience.com  
Customer Service:  
888.999.1371  
858.642.2058  
info@eBioscience.com  
Fax:  
858.642.2046

**Austria**

Technical Support:  
tech@eBioscience.com  
Customer Service:  
+43 1 796 40 40 305  
Austria@eBioscience.com  
Fax:  
+43 1 796 40 40 400

**Belgium**

Technical Support:  
tech@eBioscience.com  
Customer Service:  
+43 1 796 40 40 308  
Belgium@eBioscience.com  
Fax:  
+43 1 796 40 40 400

**France**

Technical Support:  
tech@eBioscience.com  
Customer Service:  
0 800 800 417  
France@eBioscience.com  
Fax:  
0 800 800 418

**Germany**

Technical Support:  
tech@eBioscience.com  
Customer Service:  
069 33 29 64 56  
Germany@eBioscience.com  
Fax:  
069 255 77 335

**Ireland**

Technical Support:  
tech@eBioscience.com  
Customer Service:  
+44 208 951 4482  
Ireland@eBioscience.com  
Fax:  
+44 208 952 7527

**Netherlands**

Technical Support:  
tech@eBioscience.com  
Customer Service:  
+43 1 796 40 40 308  
Netherlands@eBioscience.com  
Fax:  
+43 1 796 40 40 400

**Switzerland**

Technical Support:  
tech@eBioscience.com  
Customer Service:  
+41 21 510 1214  
Switzerland@eBioscience.com  
Fax:  
+41 21 510 1216

**United Kingdom**

Technical Support:  
tech@eBioscience.com  
Customer Service:  
+44 208 951 4482  
UK@eBioscience.com  
Fax:  
+44 207 900 1559

Customers in countries where direct sales are not available may contact their eBioscience distributor listed at [www.eBioscience.com/distributors](http://www.eBioscience.com/distributors)

[www.eBioscience.com](http://www.eBioscience.com)