

Maximizing sample quality for cryo-EM single particle analysis

Sample quality is essential for optimized, high-quality cryo-electron microscopy (cryo-EM) data for single particle analysis (SPA) and reconstruction. The sample preparation process consists of multiple variables that can be adjusted to improve specimen quality. This white paper highlights a number of available options for altering protein behavior during vitrification, which can greatly improve the outcomes of cryo-EM experiments.

Over the past 10 years, cryo-EM imaging has progressively enabled ever higher resolution structures, leading to a number of revolutionary insights. The necessary high-quality images for these results, however, rely heavily on the quality of the samples. For single particle analysis, this consists of a biological macromolecule (normally a protein) that is well-folded, homogeneous, non-aggregated, monodisperse, and preserved in vitreous ice at good density. The resulting sample should have a high number of particles in each image that are not overlapping or otherwise aggregated. If biochemical information is available about native protein structure and oligomeric state, the sample should also conform to this data.

Sample preparation is a critical part of the cryo-EM workflow, and it is therefore important to have a good understanding of each step in the process, along with selection of appropriate reagents, tools, and instruments. A typical sample preparation workflow contains the following steps:

Table 1. Vitrification workflow

Stage	Summary	Techniques
 1. Protein expression	The protein expression system can change protein abundance, as well as solubility, aggregation state, and other factors.	Use the Thermo Scientific™ Gibco™ Protein Expression System Tool
 2. Protein purification	Biochemical purification and buffer condition screening (as well as options such as cross-linking) produce homogeneous, well-folded, and monodisperse protein samples.	A range of chromatography and other tools, as seen on the Thermo Fisher Scientific protein sample preparation page including Thermo Scientific™ VitroEase™ Buffer Screening Kit
 3. Quality assessment	Assessment of purification quality, through biochemical assays and/or negative-stain TEM, helps ensure the sample is suitable for vitrification.	Stability and size assays, as well as room-temperature TEM
 4. Vitrification	The vitrification process includes many variables that can affect the density, state, and dispersion of the protein on the grid.	Vitrification device configuration along with grid selection and modification
 5. Grid screening	Identifying the highest quality grids for high-resolution data collection requires careful inspection of the grids.	Cryo-TEM



Protein purification and biochemistry

Most successful imaging protocols utilize a protein solution that is approximately 1 mg/mL of homogeneous, well-folded, monodisperse protein. If a sample does not conform naturally to this challenging set of requirements, changes to the protein purification conditions and protocol should be explored.

Adjusting the conditions in which proteins are purified, stored, and utilized is one of the fastest and simplest ways to affect the behavior and condition of the sample.^{10,21,32} Variables include (but are not limited to) pH and buffer identity, as well as salt, ligands, detergents, and small molecule additives. Buffer has been shown to improve cryo-EM samples in a number of experiments, including in the analysis of amyloid fibrils²⁹ and TMEM16.⁸ There are many commercially available screens to identify the conditions that are optimal for cryo-EM experiments, including the Thermo Scientific [VitroEase Buffer Screening Kit](#).

Table 2. Common changes to buffer components that can improve protein quality

Issue	Buffer components to vary	Potential results
Aggregated protein and/or compositional heterogeneity	Salt concentration and pH. For membrane proteins, detergents, amphipols or nanodiscs can be used	Cleaner micrographs, fewer clumps, and more usable particles per hole, as well as better 2D class averages ¹²
Air-water interface damage and/or preferred orientation	Detergent variation, particularly CHAPSO	More orientations, as well as improved 3D reconstruction isotropy ^{5,10}
Need for stabilizing additives	Small molecule additives	Improved complex stability/monodispersity, as well as higher achievable resolution ³
Unexpected interactions between the buffer and sample support	Any, but with a particular focus on detergents (when present) and small molecule additives	More consistent ice along with more even particle distribution ²⁹

Simple modifications to the protein surface, such as PEGylation,³³ will change surface charge, and can therefore change vitrification behavior. As a result, this should also be considered as an option with samples that show a tendency to aggregate.

Cross-linking can stabilize complexes for cryo-EM

Cross-linking is often used in single particle analysis when a macromolecular assembly is fragile, dynamic, or held together by weak interfaces that tend to fall apart (or rearrange) during purification, blotting, and vitrification. By introducing low, carefully titrated covalent stabilization, subunit dissociation and “breathing” motions can be reduced enough to improve particle integrity, 2D class quality, and the fraction of usable particles, especially for multi-protein/RNP complexes that otherwise show compositional heterogeneity. A key tradeoff is that cross-linking can also lock in non-native states, reduce the conformational diversity you actually want to resolve, and/or introduce subtle structural distortion, so conditions are typically screened gently and validated with activity/binding assays and negative-stain EM.

Table 3. Common cross-linking stabilization methods

Method (common name)	Summary	Typical chemistry and reagents	When to use	Common problems
GraFix (gradient fixation) ¹⁴	Rate-zonal centrifugation through a density gradient that contains a fixative	Glutaraldehyde (or other mild fixatives)	For very fragile, multi-subunit complexes; reduces dissociation while purifying	Over-fixation can blur high-resolution features; gradient components must be removed before freezing
Gentle solution cross-linking (batch crosslinking) ¹³	Adds low-dose crosslinker directly in solution, then quenches	Glutaraldehyde (or formaldehyde-like fixatives)	For quick stabilization when gradients aren't needed	Can induce aggregation when higher sample concentrations are used during reaction
AgarFix (agarose fixation) ¹	Embeds sample in an agarose drop; the crosslinker is then diffused and the sample recovered	Mild chemical crosslinkers, which are diffused into the matrix	When batch crosslinking causes aggregation; for "sticky" complexes	The recovery step adds handling time, and diffusion time/conditions need tuning
Amine-reactive NHS-ester crosslinking ²⁴	Crosslinks via lysines/N-termini using NHS esters, then quenches with Tris/glycine	BS3, DSS, or DSG (families widely used in structural biology)	For capturing weak interfaces and stabilizing assemblies; also common in integrative workflows	Can preferentially modify functional lysines, and spacer length/solubility can have an impact; additionally, hydrolysis competes with the crosslinking
Zero-length crosslinking ²⁴	Creates "no spacer" links between proximate functional groups	EDC (carbodiimide, which couples carboxyls to amines)	For mapping very close contacts; sometimes used to rigidify interfaces	Highly condition-dependent due to pH/side reactions; can also be harsh for delicate complexes

If these methods are not enough to keep the protein in a suitable state, changes to expression construct, vector, and protocol may be required in order to obtain a more soluble, well-folded macromolecule. Thermo Fisher Scientific provides extensive support for [relevant expression and purification workflows](#).



Protein quality assessment

Before vitrification and high-resolution imaging in a high-end cryo-transmission electron microscope (cryo-TEM) is attempted, it is important to establish that the expression and purification protocol has produced a sample of appropriate stoichiometry and complex stability. Ideally, this means that it consists of well-folded macromolecules in a suitably concentrated, monodisperse, and homogeneous solution. There are a wide variety of biophysical assays that can be used to assess protein state, including:

- **Differential scanning fluorometry**,^{18, 19} also known as ThermoFluor, which measures thermal unfolding (melting temperature) through the binding of an extrinsic hydrophobic dye (e.g., SYPRO Orange). These dyes often interact with detergents, so this method is not recommended for membrane proteins.
- **ProteoPlex**,⁴ which measures thermal unfolding and aggregation behavior for macromolecular complexes, utilizing a sparse chemical matrix.
- **CPM thermal shift assays**,² which uses a thiol-reactive reporter dye to measure thermal unfolding via the exposure of cysteines. This method works well with membrane proteins.
- **SEC-MALS**,²⁶ which measures absolute mass and polydispersity across SEC peaks, thereby reporting oligomeric states, as well as the presence of aggregates or polydispersity.
- **Mass photometry**,²⁷ which reports on single-molecule mass distribution via light scattering during transient surface binding. It works at low concentrations and is label-free, aiding in the detection of weak/transient assemblies, or minor populations.
- **Negative-stain electron microscopy**,²⁸ which provides direct visualization of particle integrity, morphology, and flexibility, making it one of the preferred final screening tools prior to cryo-EM. However, its low resolution, as well as its occasional staining and surface artifacts, make it unreliable for stoichiometry. Thermo Scientific™ Talos™ 12/F200C TEMs are an ideal choice for negative stain pre-screening prior to vitrification, though desktop S(T)EM devices, such as the Thermo Scientific™ Phenom Pharos™ Desktop SEM, can provide a useful alternative. Negative stained samples can be prepared with the ready-to-use Thermo Scientific™ [VitroEase Methylamine Vanadate Negative Stain](#) or the [VitroEase Methylamine Tungstate Negative Stain](#).
- **Native mass spectrometry (MS)** ³⁴ measures the mass of proteins and complexes in native-like conditions, enabling the determination of subunit composition, stoichiometry, ligand binding, and heterogeneity. It is best used for soluble complexes, protein-ligand binding, and samples where high-accuracy stoichiometry/composition is important.

The choice of assay is dependent on protein type and the characteristic being investigated (Table 4).

Table 4. Suggested assays by protein type

Protein type	Stability	Size/monodispersity	Pre-cryo-EM integrity check
Soluble globular proteins	ThermoFluor (DSF)	SEC-MALS, mass photometry	Negative stain
Multi-subunit soluble complexes	ProteoPlex	SEC-MALS, mass photometry, negative stain, native MS for stoichiometry	Negative stain
Membrane proteins	CPM	SEC-MALS	Negative stain
Weak/transient assemblies	DSF or CPM (with ligands)	Mass photometry (primary), native MS for stoichiometry/ligand-binding	Negative stain

Vitrification

In the vitrification step, the protein solution is dispensed onto a sample support, and then rapidly frozen with the aid of a plunge-freezing device. For high-resolution studies, the resulting vitrified specimen should have well-dispersed protein particles with as high a density as possible without crowding and aggregation. It should be frozen in a flat, continuous ice layer, no thicker than ~100 nm (in order to allow transmission of the electron beam). Note that variations in both the sample support and vitrification settings can have a significant impact on ice thickness, sample density, and dispersion.

Plunge freezing device

The Thermo Scientific™ [Vitrobot™ Mark IV System](#) is an established tool that offers semi-automated vitrification for simplified and routine cryo-EM sample preparation. It performs the cryo-fixation process at constant physical and mechanical conditions (i.e., temperature, relative humidity, blotting conditions, and freezing velocity). Improved stability for the cryo-fixation environment can help improve reproducibility and quality of sample preparation prior to cryo-TEM observation. A range of variables can be adjusted in the Vitrobot System in order to alter sample particle dispersion and density. The Thermo Scientific™ [VitroEase Cryo-EM Training Kit](#), which offers a detailed visual manual, is also available to help new users learn the intricacies of the vitrification process. Some default settings that work well for many proteins are given below, together with potential adjustments that can be made for specific issues (Table 5).

Table 5. Default Vitrobot settings and potential adjustments

Variable	Suggested setting	Issue	Variable to adjust
Temperature	4°C	Ice too thick	Increase blot time or force, increase temperature, reduce sample concentration
Humidity	100%	Ice too thin/empty holes	Decrease blot time or force, increase humidity
Wait time	5–10 s	Particles denatured/at air-water interface	Reduce wait or blot time, lower temperature
Blot force	0	Preferred orientation	Reduce blot time, lower temperature, Increase ice thickness, adjust blot force asymmetrically (offset)
Blot time	3–4 s	Aggregation	Lower temperature, reduce wait time or blot force
Blot count	1	Ice gradient across grid	Adjust blot offset, check pad condition, avoid excessive blot force
Blot offset	0 mm		
Drain time	0 s		
Plunge speed	Default		
Cryogen	Fresh liquid ethane		

Specimen support selection

For most cryo-EM experiments, the protein sample is dispensed onto a 3-mm-diameter grid, which is coated with a sample support film (Figure 1).






Figure 1. Diagram of a cryo-EM grid, highlighting the mesh, support film, and other layers. Images courtesy of Quantifoil Micro Tools GmbH.

Both the grid and the support film can be made from a variety of materials and with a range of different patterns. Grid selection is dependent on both the nature of the protein sample and the type of experiment. For cryo-EM, copper and gold grids are normally used. They are both conductive and radiation hard, which helps to ensure that electrons are rapidly conducted away from the sample, reducing radiation damage. 200–400 size square meshes are usually chosen, which allows for adequate sample support while still maintaining large open imaging areas. There is a broad range of options available for the choice of sample support film (Figure 2). For the highest resolution data, gold foils on gold grids are generally preferred, as these produce flatter, thinner, more consistent ice and reduce specimen movement on exposure to the beam, enabling higher resolution reconstructions.^{17, 23}




Continuous films

- Formvar/polymer-resilient, not radiation hard (<60kV scopes only)
- Carbon—more radiation hard but much more fragile



Holey/lacey carbon

- Random arrangement of random hole sizes
- Can be made in-house or bought from suppliers
- Don't work well with automation
- Good with variable-size samples



Defined hole films⁷

- Carbon (e.g. Quantifoil, C-Flat), SiO₂, SiN, polymer
- Holes of pre-defined size and arrangement
- Good with automation
- Radiation stable, but mechanically fragile
- Made with lithography
- Cryo-EM standard choice





Image courtesy of Calibre Scientific Ltd

Metal foils^{11,17,23}

Foils made by lithography from gold, palladium, amorphous nickel-titanium alloy (ANTA) and others e.g. UltrAuFoil, HexAuFoil, Au-Flat, ANTCryo

- Overcome issues with foil stability
- Produce flatter, thinner ice
- Reduce specimen movement



Silicon Nitride/SiO₂

- Example: CryoSilico
- Resilient films, very biocompatible
- Stable at higher temperatures
- Can provide modifiable surface for chemistry




Image courtesy of EM Resolutions Ltd.

Additional layers

- Graphene, GO, amorphous carbon
- Improve distribution and density by providing a surface across holes to absorb particles/provide additional support
- Can be modified to introduce additional grid functionality

Figure 2. Common options for grid support films. Unless otherwise stated, images courtesy of Quantifoil Micro Tools GmbH.

Support film hole size can also change the behavior of the grid; as a result, support films are available with a large range of hole sizes, shapes, and geometries. The universal naming convention for grids is R x/y, where R/S is the hole shape ("R" is round and "S" is square), x is the hole diameter, and y is the space between the holes. Some examples of common film geometries for single particle analysis data collection are shown in Figure 3.




Figure 3. Common support film geometries used for cryo-EM single particle analysis. Images courtesy of Quantifoil Micro Tools GmbH.

While carbon and gold supports can have holes as large as 17 μm , these do not generally support thin layers of vitreous ice. Instead, for most cryo-EM imaging, films with holes $\leq 2 \mu\text{m}$ are used. In general, smaller holes produce thinner, flatter ice, as they provide greater support. There is evidence that holes as small as 0.3 μm can be beneficial for high resolution data collection.¹⁷ However, considerations such as the number of particles per hole, and the use of multi-shot data collection protocols, mean that the most common foil geometry is R 1.2/1.3. This provides a balance between thin ice and strong support for high-resolution imaging, together with reasonable particle numbers for processing. Larger proteins and viruses may benefit from the use of R 2/2 or R 2/1 (i.e., with 2 μm holes), which improves particle numbers per hole.

The latest innovation in specimen supports for cryo-EM optimize grid and support film combinations to maximize performance. For example, HexAuFoil sample supports utilize gold foils on a gold grid; these minimize beam-induced motion in the plane of the grid, with small, 0.3 μm hole sizes that reduce strain on vitrification and minimize ice thickness, along with a hexagonal mesh grid that provides maximal sample support without compromising the open viewing area.¹⁷

Grid preparation

Both carbon and metal foils are hydrophobic, and protein solutions do not naturally spread across the grid (a process called “grid wetting”). As a result, the grids must be placed in a glow-discharger or plasma-cleaner for a few seconds (usually 10–60 seconds, depending on the model and settings) before they can be used. Longer exposure will increase hydrophilicity and grid wetting, but prolonged exposure can damage the support foil and can also lead to dry holes as the protein solution spreads too quickly and becomes too thin. Plasma cleaners can use different gases, such as argon, for this process; this can change wetting behavior and increase reproducibility compared to glow discharging in air.²¹ In addition, volatile additives, such as poly-lysine, can be added to the chamber; this modifies the film surface and the behavior of the protein, both on the grid and under the beam.



Figure 4. Illustration showing the effect of glow discharging on the behavior of a sample solution on a cryo-EM grid. The increase in surface hydrophilicity facilitates the distribution of the aqueous sample across the grid.

Grid modifications

In addition to glow discharging of the grid to improve wettability, there are a range of chemical and physical modifications that can be made to the grid surface to alter its interaction with the specimen solution, potentially improving sample stability, density, or dispersion. These include:

- **Additional ultrathin layers of continuous amorphous carbon**, which can be purchased already applied to the grids or added in-house.²¹ Many proteins are adsorbed by carbon, which prevents them from appearing within the grid holes. A continuous carbon layer allows proteins to adsorb across the holes, and simultaneously sequesters the protein away from the air-water interface, preventing denaturation. This does increase background noise, so it is only suitable for larger proteins and complexes. Additionally, this modification can have variable effects on preferred orientation, either increasing or decreasing it, depending on the sample.
- **Monolayer graphene**,⁹ usually CVD graphene, works similarly to an amorphous carbon layer, but has the advantage of improved electron transparency, so it can be used with smaller proteins. However, graphene can hold significant charge in the beam which can cause damage to the protein being imaged.
- **Graphene oxide (GO) support layers**²² are similar to graphene monolayers; however, GO flakes are stable and readily purchased from many suppliers, making preparation more straightforward. Conversely, GO flakes are rarely monolayers, or even consistent thickness, and therefore, coverage can vary.
- **Functionalized graphene/GO layers**¹⁵ have been chemically modified to introduce a defined surface charge or capture a specific protein tag. They can alter electrostatics to reduce preferred orientation, sequester particles away from the air-water interface, or reduce non-specific protein binding.
- **PEGylation/passivation of grids**¹⁶ involves the addition of PEG or PEG derivatives to the grid surface, making it more hydrophilic and “non-fouling.” It reduces non-specific adsorption onto the grids, thereby increasing partitioning of particles into the holes. Combined with PEG spacers, the technique can also offer control of sample particle orientation and/or distance from grid.
- **Affinity grids** have been chemically modified to include functional layers that bind tagged samples, sequestering protein away from the air-water interface and increasing the concentration of scarce samples. Many different varieties exist, including those with Ni-NTA²⁵ and other common tags, such as SpyCatcher.³⁰ An alternative to chemical modification is the growth of a streptavidin monolayer crystal on the surface of the grid.⁶ This is useful for biotinylated proteins, and the streptavidin background can also be easily subtracted from images in Fourier space.³¹

Table 6 provides a guide to which of the wide range of possibilities for grid modifications are likely to be the most use for a given imaging issue.

Table 6. Common sample issues and possible grid modifications to address them

Issue	Relevant grid modification	Outcome
Low particle density in holes	Ultrathin carbon, graphene/GO	Provides a protein adsorption surface across holes, increasing partitioning to holes
Air-water interface denaturation/particles are all at surface	Graphene/GO, affinity grids (Ni-NTA, streptavidin), PEG passivation	Moves particles to grid surface, away from air-water interface
Strong preferred orientation	Graphene/GO (including functionalized/charged graphene), ultrathin carbon, PEG-SAMs, affinity capture with altered tether geometry	Changes the dominant adsorption interface and/or electrostatics
Beam-induced motion is limiting resolution	Graphene/GO or change to gold film and reduce hole size	Stabilizes ice and reduces charging/motion
Scarce sample that needs concentrating	Ni-NTA lipid monolayer, streptavidin grids, GO (which often allows for lower concentrations)	Increases on-grid concentration and/or capture efficiency

Screening

With so many possible variables, a robust screening protocol is essential for optimizing the vitrification protocol. This can be done by performing diagnostic cryo-EM imaging before high-resolution data acquisition begins. This step qualitatively assesses if the sample is a promising target for 2D class average analysis while also, simultaneously, obtaining an initial low-resolution reconstruction. During this step, the sample is pre-screened to evaluate the following properties:

- Protein concentration, stability, and distribution
- Ice quality, thickness, and uniformity across the grid

Grid screening is often repetitive, testing multiple sample and vitrification options in order to identify optimal conditions.

Thermo Scientific Smart EPU Software offers fully automated grid screening to simplify this process and increase throughput.³⁵ During the screening phase, throughput can be further increased with EPU Multigrad Software, which enables the imaging of multiple grids in a single, unattended session.

If the sample is promising, a larger set of images may be acquired to facilitate further 2D and 3D analyses. At this stage, only a moderate resolution 3D map (>3 Å) is required. Thermo Scientific™ **Glacios™ Cryo-TEMs** and Thermo Scientific™ **Tundra™ Cryo-TEMs** are well suited for this screening, and offer robust, contamination-free connectivity with the higher resolution Thermo Scientific™ **Krios™ Cryo-TEM**.

References

1. Adamus K, *et al.* **AgarFix: Simple and accessible stabilization of challenging single-particle cryo-EM specimens through crosslinking in a matrix of agar.** *J Struct Biol* 207:3 (2019). doi: [10.1016/j.jsb.2019.07.004](https://doi.org/10.1016/j.jsb.2019.07.004)
2. Alexandrov AI, *et al.* **Microscale fluorescent thermal stability assay for membrane proteins.** *Structure* 16:3 (2008). doi: [10.1016/j.str.2008.02.004](https://doi.org/10.1016/j.str.2008.02.004)
3. Basanta B, *et al.* **A case for glycerol as an acceptable additive for single-particle cryoEM samples.** *Acta Cryst D* 78:1 (2022). doi: [10.1107/s2059798321012110](https://doi.org/10.1107/s2059798321012110)
4. Chari A, *et al.* **ProteoPlex: stability optimization of macromolecular complexes by sparse-matrix screening of chemical space.** *Nat Methods* 12 (2015). doi: [10.1038/nmeth.3493](https://doi.org/10.1038/nmeth.3493)
5. Chen J, *et al.* **Eliminating effects of particle adsorption to the air/water interface in single-particle cryo-electron microscopy: Bacterial RNA polymerase and CHAPSO.** *J Struct Biol X* 1 (2019). doi: [10.1016/j.yjsbx.2019.100005](https://doi.org/10.1016/j.yjsbx.2019.100005)
6. Cookis T, *et al.* **Streptavidin-Affinity Grid Fabrication for Cryo-Electron Microscopy Sample Preparation.** *J Vis Exp* 202:10.3791/66197 (2023). doi: [10.3791/66197](https://doi.org/10.3791/66197)
7. Ermantraut E, Wohlfart K, and Tichelaar W. **Perforated support foils with pre-defined hole size, shape and arrangement.** *Ultramicroscopy* 74:1-2 (1998). doi: [10.1016/S0304-3991\(98\)00025-4](https://doi.org/10.1016/S0304-3991(98)00025-4)
8. Feng S, *et al.* **Cryo-EM Studies of TMEM16F Calcium-activated ion channel suggest features important for lipid scrambling.** *Cell Rep* 28:2 (2019). doi: [10.1016/j.celrep.2019.06.023](https://doi.org/10.1016/j.celrep.2019.06.023)
9. Grassetti AV, May MB, and Davis JH. **Application of monolayer graphene to cryo-electron microscopy grids for high-resolution structure determination.** *J Vis Exp* 201:e66023 (2023). doi: [10.3791/66023](https://doi.org/10.3791/66023)
10. Han BG, *et al.* **Challenges in making ideal cryo-EM samples.** *Curr Opin Struct Biol* 81 (2023). doi: [10.1016/j.sbi.2023.102646](https://doi.org/10.1016/j.sbi.2023.102646)

11. Huang X, *et al.* Amorphous nickel titanium alloy film: A new choice for cryo electron microscopy sample preparation. *Prog Biophys Mol Biol* 156 (2020). doi: [10.1016/j.pbiomolbio.2020.07.009](https://doi.org/10.1016/j.pbiomolbio.2020.07.009)
12. Kampjut D, Steiner J, and Sazanov LA. Cryo-EM grid optimization for membrane proteins. *iScience* 24:3 (2021). doi: [10.1016/j.isci.2021.102139](https://doi.org/10.1016/j.isci.2021.102139)
13. Kang JY, *et al.* Structural basis of transcription arrest by coliphage HK022 Nun in an Escherichia coli RNA polymerase elongation complex. *eLife* 6:e25478 (2017). doi: [10.7554/eLife.25478](https://doi.org/10.7554/eLife.25478)
14. Kastner B, *et al.* GraFix: sample preparation for single-particle electron cryomicroscopy. *Nat Methods* 5 (2008). doi: [10.1038/nmeth1139](https://doi.org/10.1038/nmeth1139)
15. Liu N and Wang HW. Graphene in cryo-EM specimen optimization. *Curr Opin Struct Biol* 86 (2024). doi: [10.1016/j.sbi.2024.102823](https://doi.org/10.1016/j.sbi.2024.102823)
16. Meyerson JR, *et al.* Self-assembled monolayers improve protein distribution on holey carbon cryo-EM supports. *Sci Rep* 4:7084 (2014). doi: [10.1038/srep07084](https://doi.org/10.1038/srep07084)
17. Naydenova K, Jia P, and Russo CJ. Cryo-EM with sub-1 Å specimen movement. *Science* 370:6513 (2020). doi: [10.1126/science.abb7927](https://doi.org/10.1126/science.abb7927)
18. Niesen FH, Berglund H, and Vedadi M. The use of differential scanning fluorimetry to detect ligand interactions that promote protein stability. *Nat Protoc* 2 (2007). doi: [10.1038/nprot.2007.321](https://doi.org/10.1038/nprot.2007.321)
19. Pantoliano MW, *et al.* High-density miniaturized thermal shift assays as a general strategy for drug discovery. *J Biomol Screen* 6:6 (2001). doi: [10.1177/108705710100600609](https://doi.org/10.1177/108705710100600609)
20. Michael W, *et al.* High-Density Miniaturized Thermal Shift Assays as a General Strategy for Drug Discovery. *J Biomol Screen* 6:6 (2001). doi: [10.1177/108705710100600609](https://doi.org/10.1177/108705710100600609)
21. Passmore LA and Russo CJ. Chapter Three - Specimen Preparation for High-Resolution Cryo-EM. *Meth Enzym* 579 (2016). doi: [10.1016/bs.mie.2016.04.011](https://doi.org/10.1016/bs.mie.2016.04.011)
22. Palovcak E, *et al.* A simple and robust procedure for preparing graphene-oxide cryo-EM grids. *J Struct Biol* 204:1 (2018). doi: [10.1016/j.jsb.2018.07.007](https://doi.org/10.1016/j.jsb.2018.07.007)
23. Russo CJ and Passmore LA. Ultrastable gold substrates for electron cryomicroscopy. *Science* 346:6215 (2014). doi: [10.1126/science.1259530](https://doi.org/10.1126/science.1259530)
24. Schmidt C and Urlaub H. Combining cryo-electron microscopy (cryo-EM) and cross-linking mass spectrometry (CX-MS) for structural elucidation of large protein assemblies. *Curr Opin Struct Biol* 46 (2017). doi: [10.1016/j.sbi.2017.10.005](https://doi.org/10.1016/j.sbi.2017.10.005)
25. Skrajna A, *et al.* Nickel-NTA lipid-monolayer affinity grids allow for high-resolution structure determination by cryo-EM. *J Struct Biol* 217:4 (2025). doi: [10.1016/j.jsb.2025.108253](https://doi.org/10.1016/j.jsb.2025.108253)
26. Some D, *et al.* Characterization of Proteins by Size-Exclusion Chromatography Coupled to Multi-Angle Light Scattering (SEC-MALS). *J Vis Exp* 148:e59615 (2019). doi: [10.3791/59615](https://doi.org/10.3791/59615)
27. Sonn-Segev A, *et al.* Quantifying the heterogeneity of macromolecular machines by mass photometry. *Nat Commun* 11:1772 (2020). doi: [10.1038/s41467-020-15642-w](https://doi.org/10.1038/s41467-020-15642-w)
28. Stark H and Chari A. Sample preparation of biological macromolecular assemblies for the determination of high-resolution structures by cryo-electron microscopy. *Microscopy* 65:1 (2016). doi: [10.1093/jmicro/dfv367](https://doi.org/10.1093/jmicro/dfv367)
29. Valli D, *et al.* Improving cryo-EM grids for amyloid fibrils using interface-active solutions and spectator proteins. *Biophys J* 123:6 (2024). doi: [10.1016/j.bpj.2024.02.009](https://doi.org/10.1016/j.bpj.2024.02.009)
30. Wang F, *et al.* General and robust covalently linked graphene oxide affinity grids for high-resolution cryo-EM. *Proc Natl Acad Sci USA* 117:39 (2020). doi: [10.1073/pnas.2009707117](https://doi.org/10.1073/pnas.2009707117)
31. Wang L and Sigworth FJ. Liposomes on a streptavidin crystal: a system to study membrane proteins by cryo-EM. *Methods Enzymol* 481 (2010). doi: [10.1016/s0076-6879\(10\)81007-9](https://doi.org/10.1016/s0076-6879(10)81007-9)
32. Wang L and Zimanyi CM. Cryo-EM sample preparation for high-resolution structure studies. *Acta Cryst F* 80 (2024). doi: [10.1107/s2053230x24002553](https://doi.org/10.1107/s2053230x24002553)
33. Zhang Z, *et al.* Improving particle quality in cryo-EM analysis using a PEGylation method. *Structure* 29:10 (2021). doi: [10.1016/j.str.2021.05.004](https://doi.org/10.1016/j.str.2021.05.004)
34. Olinares PDB *et al.* Native Mass Spectrometry-Based Screening for Optimal Sample Preparation in Single-Particle Cryo-EM. *Structure* 29:2 (2021). doi: [10.1016/j.str.200.11.001](https://doi.org/10.1016/j.str.200.11.001)
35. Grollios F *et al.* Fully Automated Screening Workflow with Smart EPU. *Microscopy and Microanalysis* 31:S1 (2025). doi: [10.1093/mam/ozaf048.527](https://doi.org/10.1093/mam/ozaf048.527)

 Learn more about sample preparation for cryo-electron microscopy at thermofisher.com/cryoem