remel**BactiDrop™ India Ink EN**

REF R21518

**INTENDED USE**

BactiDrop India Ink is a stain for the microscopic examination of encapsulated organisms in clinical samples. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having fungal infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

SUMMARY AND EXPLANATION

Organisms that possess a polysaccharide capsule exhibit a halo around the cell against a black background. India ink has been reported useful in observing *C. neoformans* in cerebrospinal fluid, which aids in the diagnosis of meningitis.¹

PRINCIPLE

In the India ink method, the polysaccharide capsule displaces the colloidal carbon particles of the ink and appears as a clear halo around the microorganism.

REAGENT

Pelikan™ Drawing Ink #17 Black

PRECAUTIONS

WARNING! May impair fertility. May cause harm to the unborn child. May cause eye, skin, and respiratory tract irritation.

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully. Refer to Safety Data Sheet on company website for detailed information on reagent chemicals.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction, do not use the device

STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 20-25°C until used. Do not freeze or overheat.

PRODUCT DETERIORATION

This product should not be used if (1) there is evidence of dehydration, (2) the color has changed, (3) the expiration date has passed, or (4) there are other signs of deterioration. The expiration date applies to the product in its intact container when stored as directed.

SPECIMEN COLLECTION, STORAGE, TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{2,3}

MATERIALS PROVIDED

50 BactiDrop India Ink – glass ampule enclosed in a protective plastic sheath with a dropper dispenser tip (0.75/ampule)

1 Ampule crusher (reusable)

MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Glass slides, coverslips, (7) Microscope, (8) BactiDrop Potassium Hydroxide (10%) (R21524).

PROCEDURE

Place dropper in the assembled, reusable ampule crusher provided. Hold the dropper/crusher in an upright position and lightly tap the bottom to dislodge any bubbles that may have formed. Grasp the middle of the dropper/crusher with the thumb and forefinger, and with the tip pointing away, press gently to crush the ampule. Invert dropper and squeeze slightly to dispense in a dropwise fashion. Discard remaining portion of partially used ampule at end of workday.

Test Procedure:

1. Specimens (sputum or pus) may be digested in a potassium hydroxide mount, which eliminates most cells and other artifacts that may be misread as *Cryptococcus*.
2. Mix the specimen (pus, exudate, tissue, sputum, or sediment of centrifuged spinal fluid) with a small drop of India Ink on a clean slide. Cover with a coverslip. The mount should be thin and a hazy brown color, not thick and black.
3. Scan under low power, with reduced lighting. Switch to high power (400 x) to examine for the presence of encapsulated cells.

INTERPRETATION

C. neoformans appears as a spherical, oval or elongated, single or multiple budding, thick-walled, yeast-like organism 2-15µm in diameter, usually surrounded by a wide, refractile polysaccharide capsule. The capsule appears as a clear halo that surrounds the yeast cell or lies between the cell wall

and the surrounding black mass of India ink particles.^{4,5}

QUALITY CONTROL

All lot numbers of BactiDrop India Ink have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL

Cryptococcus gattii
ATCC® 34877

Candida albicans
ATCC®10231

RESULTS

Clear halo observed
against a dark
background

No halo observed

LIMITATIONS

1. India Ink preparations can be diagnostic; however, they are presumptive and may lack sensitivity. Further testing, such as cryptococcal latex agglutination and culture, is required for definitive identification.^{4,5}
2. Red and white blood cells, glial cells, macrophages and other artifacts may resemble cryptococcal yeasts in an India ink preparation. Further testing is required for identification.^{4,5}
3. Some strains of *Cryptococcus neoformans* may not possess detectable capsules making culture essential.⁴

BIBLIOGRAPHY

1. Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy. 1991. Manual of Clinical Microbiology. 5th ed. ASM, Washington, D.C.
2. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, and D.W. Warnock. 2011. Manual of Clinical Microbiology. 10th ed. ASM Press, Washington, D.C.
3. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.
4. Baron, E.J. and S.M. Tenegold. 1990. Bailey and Scott's Diagnostic Microbiology. 8th ed. The C.V. Mosby Co., St. Louis, MO.
5. Rippon, J.W. 1982. Medical Mycology, The Pathogenic Fungi and The Pathogenic Actinomycetes. 2nd ed. W.B. Saunders Co., Philadelphia, PA.

Symbol Legend

	Catalog Number
	In Vitro Diagnostic Medical Device
	Batch Code (Lot Number)
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
	Use By (Expiration Date)
	Contains sufficient for <n> tests
	Do not use if packaging damaged
	Authorized European Representative
	Manufacturer
	UK Conformity Assessment
	European Conformity Assessment

Remel Inc.
12076 Santa Fe Trail Drive
Lenexa, KS 66215, USA
www.thermofisher.com/microbiology
(800) 255-6730
International: (913) 888-0939

For technical information contact your local distributor.

©2022 Thermo Fisher Scientific Inc. All rights reserved
BactiDrop is a trademark of Thermo Fisher Scientific and its subsidiaries.

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

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Pelikan is a trademark of Pelikan Vertriebsgesellschaft MBH & CO.KG
Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, EE. UU.