

Specifications

Product Description : Sodium Acetate Anhydrous SQ 500G
Product Code : Q20125

Sr.No	Description	Unit	Standards
1	Description : White hygroscopic powder		Should Comply
2	Solution: 10% w/v (in water) : Almost clear & colourless		Should Comply
3	Assay	%	98 Min
4	Chloride (Cl)	%	0.01 Max
5	Heavy metals (as Pb)	%	0.002 Max
6	Iron (Fe)	%	0.002 Max
7	Potassium (K)	%	0.1 Max
8	Sulphate (SO4)	%	0.01 Max

Quality Assurance

The document has been produced electronically and is valid without signature.

Registered Office : Thermo Fisher Scientific India Pvt. Ltd. 403-404,B-wing, Delphi, Hiranandani Business Park, Powai, Mumbai-400076, INDIA.
Tel : +91-22-67162200/66803000, CIN-U73100MH2000PTC126872.