

	Document Owner Department: QC	BT-SPEC-0612
		Page 1 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
CEPHALEXIN CL30 CT0007B		

Product Details

Product Name	CEPHALEXIN
Disc Code	CL30
Disc Concentration	30µg
CT Code	CT0007B
Shelf Life (years)	3
CLSI Recommended Concentration (inc Vet Standard)	No
EUCAST Recommended Concentration	Yes

Microbiological Performance Criteria

EUCAST		
Organism	ATCC/NCTC Number	Range (mm)
<i>Escherichia coli</i>	ATCC®25922/NCTC12241	15 – 21

References

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests;

- EUCAST Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST

Other additional challenging strains may be employed.

The information given above is believed to be correct, however both the information and the product are offered without warranty for any specific use, other than that specified in the current Oxoid manual.

	Document Owner Department: QC	BT-SPEC-0612
		Page 2 of 2
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
CEPHALEXIN CL30 CT0007B		

Revision History

Section / Step	Description of Change	Reason for Change	Reference
Microbiological Performance Criteria	Removal of BSAC	Change control	BT-CC-2142
Entire document/ References	Update to current format. Update to current performance standards	Change control	BT-CC-2760