

Veterinary Antimicrobial Susceptibility Testing

Reliable results for confident animal care

Accurate antimicrobial susceptibility testing for aquatic animals

The Thermo Scientific™ Sensititre™ Veterinary Aquatic Animals AQUA2 Plate delivers accurate, minimum inhibitory concentration (MIC) results utilizing broth microdilution—the gold-standard method¹ for antimicrobial susceptibility testing (AST).

Designed by industry experts and meeting the latest CLSI guidelines, this new standard plate offers a comprehensive range of aquatic animal focused antimicrobials on a single format, ensuring accurate MIC determination to support antimicrobial resistance diagnostics, surveillance, and stewardship efforts.


The Sensititre Veterinary Aquatic Animals AQUA2 Plate is compatible to be read with the Thermo Scientific™ Sensititre™ Manual Viewbox and Thermo Scientific™ Sensititre™ Vizion™ readers, supported by the Thermo Scientific™ Sensititre™ SWIN Software System for result interpretation.

1. Gram negative anaerobe susceptibility testing in clinical isolates using Sensititre and Etest methods.
C. Hughes, C. Ashhurst-Smith, J.K. Ferguson. Pathology Volume 50, Issue 4, June 2018.



Find out more about our Veterinary testing solutions or contact us for a quote by visiting our website at [thermofisher.com/vet-microbiology](https://www.thermofisher.com/vet-microbiology)

thermo scientific



Antimicrobics

AMI	Amikacin
AMP	Ampicillin
TAZ	Ceftazidime
ENRO	Enrofloxacin
FFN	Florfenicol
GEN	Gentamicin
MERO	Meropenem
NEG	Negative Control
OXO	Oxolinic Acid
OXY	Oxytetracycline
POS	Positive Control
SXT	Trimethoprim / Sulfamethoxazole

Product description	Quantity	Product code
Consumables		
Thermo Scientific™ Sensititre™ Veterinary Aquatic Animals AQUA2 Plate	10 plates/pack	AQUA2
Thermo Scientific™ Sensititre™ Demineralized water	100 x 5 mL tubes/pack	T3339
Thermo Scientific™ Sensititre™ Mueller Hinton Broth	100 x 11 mL tubes/pack	T3462
Thermo Scientific™ Sensititre™ Mueller Hinton Broth	100 x 5 mL tubes/pack	T3462-05
Thermo Scientific™ Sensititre™ Veterinary Fastidious MHF Medium (MHF-Y)	10 x 11 mL tubes/pack	T3461
Thermo Scientific™ Sensititre™ Mueller Hinton Broth with lysed horse blood	10 x 11 mL tubes/pack	CP112-10
Instruments		
Thermo Scientific™ Sensititre™ Nephelometer	Each	V3011
Thermo Scientific™ Sensititre™ AIM™	Each	V3020
Thermo Scientific™ Sensititre™ Vizion	Each	V2021
Thermo Scientific™ Sensititre™ Manual Viewbox	Each	V4007
Thermo Scientific™ Sensititre™ Electronic Multichannel Pipette	Each	V4009
Thermo Scientific™ Sensititre™ SWIN PC Global Kit	Each	SW4000GBL*
Thermo Scientific™ Sensititre™ SWIN Software System	Each	SW4000**
Thermo Scientific™ Sensititre™ SWIN Software Epidemiology Module	Each	SW1203
Accessories		
Thermo Scientific™ 0.5 Polymer McFarland Standard	Each	E1041
Thermo Scientific™ Doseheads	100/box	E3010

*Available worldwide outside US. **Available in the US.

Find out more about our Veterinary testing solutions or contact us for a quote by visiting our website at thermofisher.com/vet-microbiology

thermo scientific