

Antimicrobial Susceptibility Testing


Discover a new MIC testing panel for non-fastidious Gram-negative organisms

The Thermo Scientific™ Sensititre™ EUMDR2F broth microdilution panel

- Following the latest EUCAST breakpoint guidelines and designed in collaboration with industry experts to support a comprehensive susceptibility testing solution.
- Minimum inhibitory concentration (MIC) determination with broth microdilution methodology, including new antimicrobials such as Aztreonam/Avibactam, Cefepime/Emnetazobactam, and Sulbactam/Durlobactam.
- Includes a broad range of relevant second line antimicrobials, including Ceftazidime/Avibactam, Colistin, Eravacycline, Imipenem/Relebactam, and Meropenem/Vaborbactam.
- For *in vitro* diagnostics use.
- Provides a consolidated design for susceptibility testing of relevant non-fastidious organisms, supporting clinical microbiology workflows.
- Integrate into fully automated, semi-automated, or manual workflows; compatible with the Thermo Scientific™ Sensititre™ ARIS HiQ™ System, OptiRead™, Vizion™ Digital MIC Viewing System, or manual read methods.
- Room-temperature storage throughout 24-month shelf life supports simplified inventory management.



Thermo Scientific™ Sensititre™ EUMDR2F



Antimicrobics

AMI	Amikacin
AZT	Aztreonam
AZA	Aztreonam / Avibactam Constant 4
FEP	Cefepime
FPE	Cefepime / Enmetazobactam Constant 8
CZA	Ceftazidime / Avibactam Constant 4
C/T	Ceftolozane / Tazobactam Constant 4
COL	Colistin
ERV	Eravacycline
IMR	Imipenem / Relebactam Constant 4
MERO	Meropenem
MEV	Meropenem / Vaborbactam Constant 8
NEG	Negative Control
P/T4	Piperacillin / Tazobactam Constant 4
POS	Positive Control
SUD	Sulbactam / Durlobactam Constant 4
TGC	Tigecycline
TOB	Tobramycin

Product description	Pack size	Product code
Thermo Scientific™ Sensititre™ EUMDR2F	10 / box	EUMDR2F
Thermo Scientific™ Sensititre™ Demineralized Water	100 x 5 mL tubes / pack	T3339
Thermo Scientific™ Sensititre™ Mueller Hinton Broth	100 x 11 mL tubes / pack	T3462
Thermo Scientific™ Culti-Loops™ Escherichia coli ATCC™ 25922™	5 loops / pack	R4607050
Thermo Scientific™ Culti-Loops™ Pseudomonas aeruginosa ATCC™ 27853™	5 loops / pack	R4607060

* Not all products are available for sale in all territories. Please inquire.

To learn more, visit thermofisher.com/sensititre

thermo scientific